Surgical perfection. For life.

## **INSTRUCTIONS FOR USE AND CARE**

**EN**

Stille Surgical and reusable Instruments,  
product group 101-119

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## INSTRUCTIONS FOR USE AND CARE

### Stille Surgical and reusable Instruments, product group 101-119

#### INTENDED USE

STILLE surgical instruments are intended as standard instruments for procedures in open surgery, and for the intended use in medical specialty. The instruments may only be used for their intended purpose by trained and qualified operators.

#### APPLICABLE FOR

These instructions are applicable for STILLE surgical instruments like: Scissors, Forceps, Artery Forceps, Needle Holders, Bone Cutters and Rongeurs, Clamps, Retractors, Curettes, Raspatories, Elevators, Dissectors, Rasps, Osteotomes, Chisels, Awls, Mallets and Micro Surgical Instruments.

#### CAUTION

It is impossible to foresee all potential risks with the product and the user shall be responsible for correct use and selection of instruments for each application. Used for a task other than what it is intended for may result in injury or a damaged instrument.

Any serious incident that occurs with the device covered by these IFU should be reported to STILLE AB and the competent authority of the Member State in which the user and/or patient is established.

STILLE AB cannot accept any liability for immediate or consequential damages caused by incorrect use or by inappropriate cleaning, sterilization and maintenance of the instruments. Use of a device for a task other than what it is intended for may damage the device and void the product warranty.

Repeated processing has minimal effect on the instruments' lifetime and there is no limitation regarding number of processes/cycles. End of life is normally determined by wear and damage due to use.

Handle the instruments carefully. This applies to instruments in general but particularly to micro instruments. Coated STILLE Diamond instruments may be sensitive for diathermy. Avoid direct contact with diathermy electrode as this may cause damage to the surface coating.

The instrument's performance characteristics is designed to be in line with its intended use. More detailed information regarding specific instruments can be requested from STILLE AB when needed.

Never immerse the instruments in a saline solution as this leads to pitting and stress corrosion.

Cleaning or sterilizing with any method containing hydrogen peroxide is not recommended for black coated instruments (e.g. SuperCut scissors).

Do not allow instruments to remain under mechanical stress since this increases the risks for corrosion and damage. Instruments with ratchets should be left open. Microsurgical instruments should never be laid on top of each other.

Only the cleaning, disinfection and sterilization processes described in these IFU have been validated.

The instruments are intended for temporary use (<60 min. under normal conditions) according to MDD 93/42/EEC,

appendix 9, rule 6 (according to MDR 2017/745 appendix 8, chapter 3, rule 6).

The marking is a vital part of a medical device. It's important to make sure that the original marking is intact and in good shape. Due to traceability and warranty purpose, the article code, lot-code, STILLE-logo and CE-mark must not be removed. This is of particular importance to make sure when servicing the instruments.

If instruments are repaired by any service supplier not authorized by STILLE AB, all warranties are void.

#### CONTRAINDICATIONS

Usage outside of the intended use. Usage of non-insulated instruments for diathermy procedures.

#### BEFORE FIRST USE OF NEW INSTRUMENTS

New instruments are delivered non-sterile. To be used in clinical settings, they must be cleaned, disinfected, inspected and sterilized before use, by following the procedure as described in this IFU. This is also valid for instruments returned from service/repair.

New instruments require particular attention when it comes to reprocessing / storage.

Staining and discoloring may occur in the vicinity of the laser marking following cleaning and disinfecting. Instruments which have been cleaned in the same cleaning and disinfecting unit and have been in circulation for a longer period do not exhibit these changes. This can be a consequence of the still very thin passivated layer on the new instruments.

The following instructions should be heeded to avoid problems of this nature:

- Always remove new instruments from the transportation packaging (closed plastic packaging can lead to the formation of condensation if exposed to temperature fluctuations).
- New instruments should not be stored in the immediate vicinity of chemicals (chemical vapors can trigger corrosion on the instruments).

Please contact STILLE if, despite several cleaning and disinfecting operations, discoloring is still evident in the vicinity of the marking.

#### POINT OF USE CARE

- Wipe off blood and/or debris from device throughout the surgical procedure, and as soon as possible after, to prevent it from drying onto the surface. Use dedicated sterile wipes.

#### PREPARATIONS FOR DECONTAMINATION

- Devices should be reprocessed as soon as it is reasonably practical following use. A maximum time of 6 hours after use is normally referred to as maximum time.
- Hinged and boxlock instruments must be cleaned in open position. Take instruments apart as much as possible. No STILLE-instruments are intended to be opened by using of a tool. The joint of bone rongeurs and cutters can be kept open by the use of STILLE Teflon Pin (Article code 119-78800).
- Soak and/or rinse heavily soiled instruments prior to cleaning to loosen any dried soil or debris. Use an enzymatic or low-alkaline (pH≤12) detergent to soak devices. Follow the detergent manufacturer's instructions for use for correct time,

temperature, water quality and concentration. Use cold tap water to rinse devices.

### MANUAL AND MECHANICAL CLEANING AND DISINFECTION

Only use manual and mechanical cleaning if there is no automated process available. Automated cleaning and disinfection are preferred.

Detergent: Enzymatic or low-alkaline (pH≤12) detergent suitable for manual cleaning. Follow producers' recommendations for concentration.

Equipment: soft bristle brushes, ultrasonic cleaner, lint-free cloths.

Water quality: At least drinking water quality if not otherwise specified.

1. Ensure that each instrument is in an open position.
2. Rinse soiled devices under running cold tap water for a minimum of 1 minute. Use a soft-bristle brush to assist in the removal of gross soil and debris. Don't use steel wool.
3. Soak device in fresh lukewarm enzymatic or low-alkaline (pH≤12) detergent solution for minimum 10 minutes. Follow the detergent manufacturer's instructions for use for correct time, temperature and concentration.
4. Rinse device with cold water for at least 30 seconds. Be sure to meticulously and vigorously flush parts that are hard to reach.
5. Prepare a fresh detergent solution using a non-foaming enzymatic or low-alkaline (pH≤12) detergent suitable for ultrasonic cleaning, with a temperature of 30–40 °C. Follow the detergent manufacturer's instructions for use for correct concentration. When using ultrasound for cleaning, instruments should be opened and completely covered with cleaning solution. Be careful to protect sharp instrument tips.
6. Clean the instruments ultrasonically for a minimum of 10 minutes, using a frequency of 35–45 kHz, 150W.
7. Prepare a new lukewarm enzymatic or low-alkaline (pH≤12) detergent solution. Use a soft bristle brush to gently scrub the devices below the water surface, until all visible dirt has been removed. Be extra careful with teeth, serrations, serrated handle parts and other parts that are more difficult to clean. Make sure to expose all areas to the detergent solution.
8. Rinse devices thoroughly with deionized (DI) or purified (PURW) water for at least 30 seconds. Be sure to meticulously and vigorously flush parts that are hard to reach. Make sure to expose all areas to the rinsing water.
9. Visually inspect the device. If any soil or chemicals remain on the device, repeat steps 2–8.
10. Perform a final rinse on the device using DI or PURW water, minimum 10 seconds.
11. Dry the device using a clean, soft, lint-free cloth or clean compressed air.

### DISINFECTION AFTER MANUAL AND MECHANICAL CLEANING:

1. Soak or apply plentifully of ethanol (70%) on device, covering all surfaces. Actuate joints, handles and other movable device features to expose all areas to the ethanol, if applicable.
2. Leave to dry on a clean lint-free cloth or similar.

### AUTOMATED CLEANING AND DISINFECTION:

Detergent: Enzymatic or low-alkaline (pH≤12) detergent suitable for cleaning in washer/disinfector. Follow producers' recommendations for concentration.

Equipment: Washer/disinfector validated according to EN ISO 15883.

Water quality: At least drinking water quality for cleaning and demineralized water for disinfection and final rinse.

Preparations for decontamination (as described above) must be performed prior to automated cleaning.

1. If any gross soil is visible on the devices, rinse device under running cold tap water for a minimum of 1 minute. Remove gross soil using a soft bristle brush. Be extra careful with teeth, serrations, serrated handle parts and other parts that are more difficult to clean.
2. Expose all areas of the instrument to the washing. Put all hinged instruments in open position, preferably on a support holding the instrument in an angled position.
3. Put the devices in an appropriate basket for the disinfector and run them through a standard cleaning cycle in the disinfector for devices.
4. The instruments should be removed from the machine immediately after finished program to prevent residual moisture from causing corrosion.

The following minimum parameters are crucial for accurate automated cleaning:

Step	Description
1. Pre rinse	2 min (x2), cold water
2. Cleaning	Minimum 10 minutes dish washing with enzymatic or mildly basic detergent at a concentration and temperature recommended by the detergent producer, however minimum 55 °C.
3. Post rinse	1 min (x2), warm water
4. Final rinse/ disinfection	Use a rinse aid at a concentration and temperature recommended by the detergent producer. Disinfection with minimum 1 min 90°C, final rinse with demineralized or pure water.
5. Drying	15 min hot air 110°C

### Drying

If a dry cycle is not included in the automated cleaning:

- Dry each device thoroughly inside and out with a clean, soft, lint-free cloth, to prevent rusting and malfunction.
- Open and close devices so that all areas are reached.

## INSPECTION

Following to the cleaning/disinfection, after cooling down to room temperature, carefully examine each instrument for cleanliness, breaks, cracks, discoloration, corrosion, deformations and malfunctions. It is especially essential to check areas such as tips, teeth, cutting edges, ratchets, joints, screws and guide pin. On insulated forceps, it is critical to examine the insulation coating for cracks and loose ends of insulation. Instruments that are worn out, corroded, deformed or damaged in any other way must be sorted out.

## LUBRICATING

After the cleaning process, lubricate surfaces by applying a thin layer of oil where metal is in contact with metal (e.g. scissors blades, box locks, screw joints). Only use a medical grade white oil intended for surgical instruments. Remove excess oil to avoid oily instruments and stains on autoclave packaging.

## FUNCTION TEST

Carefully examine and check the function of each individual instrument before packaging.

## PACKAGING

Use standardized packaging for sterilized instruments in sterile bags acc. to EN ISO 11607-1. Any divergences from the recommended packaging must be validated.

Avoid overstressing the devices, such as placing heavy objects on top of the forceps.

## STERILIZATION:

Equipment: autoclave, validated according to EN ISO 17665-1, autoclaving bags/wrapping material for steam sterilization.

## PROCEDURE:

- Wrapped instruments
- Minimum 3 minutes exposure to 134 °C steam (min 132 °C – max 137 °C)
- A minimum of 3-pulse fractionated vacuum cycle
- Drying for a minimum of 20 minutes

If using a sterilization container system, make sure to follow the manufacturer's instructions. An extended drying cycle may be needed to ensure proper drying. Remaining humidity in the container may cause discolouration or corrosion on the instruments.

## STORAGE AND TRANSPORTATION

After sterilization, devices should remain in sterilization packaging and stored in a clean, dry and moisture free area. Protect tips, sharp edges etc. with tip protectors. Make sure that no chemicals are close to or in the storage area. Use closed containers to the reprocessing area to avoid any damage on the instruments or contamination to the environment.

## MATERIAL

### Stainless steel instruments

All STILLE surgical instruments are made according to EN ISO 7153-1:2016. The high-grade stainless steel used for manufacturing STILLE instruments create a passive layer as protection of the surface. The steel however is only to a certain extent resistant against chloride ions and water containing

salts and minerals. Chloride ions mainly cause pitting and can also cause stress corrosion cracking. The greatest danger is water in which considerable quantities of common salts (sodium chloride) are dissolved.

### Titanium instruments

STILLE micro titanium instruments are made by the titanium alloy Ti6Al4V and can be handled and treated like steel instruments and no special precautions are needed. The handles on STILLE micro titanium instruments are anodised in blue colour for identification purpose.

## WARRANTY

STILLE AB warrants to the purchaser that its surgical hand instruments shall be free from defects in material and workmanship for a period of thirty (30) years after the date of STILLE AB's invoice. STILLE AB will repair or replace a damaged instrument free of charge if the damage is found to have been caused by a defect in material or workmanship. The warranty does not include normal wear and tear. The warranty is subject to the conditions that the instrument is used for the purpose intended, is cared for and stored according to our printed instructions and that any claims are made in writing within 30 days of the date the damage occurred. If instruments are repaired by any service supplier not authorized by STILLE AB, all warranties are void.

## RETURN AND DISPOSAL OF PRODUCT

If returned for service/repair, the instruments must be cleaned and disinfected according to these instructions before sending back to STILLE AB or its authorized service partner. Non-processed instruments will not be serviced. Pack the instruments safely to protect them from any damage.

For disposal of defect or out-of-date instrument, disinfect it according to these instructions. The material is recyclable as metal.

## EXPLANATION OF USED SYMBOLS

	Catalogue number
	Batch code
	Manufacturer
	Date of manufacture
	Refer to instructions for use
	non-sterile product
	CE-mark

**Rx Only (USA)** Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner

Remaining risks weighed against the benefits of the products have been treated and minimized in the risk management handling according to EN ISO 14971. New possible risks are continuously evaluated in the PMS process.

**MATERIALS AND MACHINES USED FOR  
THE VALIDATION STUDY**

Ultrasound bath: FisherBrand FB15054, 37 kHz

Detergent: Getinge Clean Enzymatic detergent 2 ml/l

Washer-disinfector: Steelco Lab600

Steam autoclave: Systec 60-DX

**APPLICABLE STANDARDS****EN ISO 17664**

Processing of healthcare products – Information to be provided by the medical device manufacturer

**EN-ISO 17665-1**

Sterilization of healthcare products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices

**EN-ISO 7153-1**

Surgical Instruments – Materials – Part 1: Metals

**EN 285**

Sterilization – Steam sterilizers – Large sterilizers

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