Surgical perfection. For life.

INSTRUCTIONS FOR USE AND CARE

EN

Stille Plaster Cast Remover,
product group 413

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SE

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INSTRUCTIONS FOR USE AND CARE

STILLE Plaster Cast Remover, product group 413

INTENDED USE

STILLE Plaster Cast Remover is to be used at the orthopaedic receptions to remove plaster cast. The Saw consists of handles, plastic profiles and wire saws. Sawing effect is obtained by moving the wire back and forth. The wire and profile are for single use only. The products may only be used for their intended purpose by trained and qualified operators.

APPLICABLE FOR

These instructions are applicable for STILLE Plaster Cast Removers only.

CAUTION

It is impossible to foresee all potential risks with the product and the user shall be responsible for correct use and selection of device. Used for a task other than what it is intended for may result in injury or a damaged device.

Any serious incident that occurs with the device covered by these IFU should be reported to STILLE AB and the competent authority of the Member State in which the user and/or patient is established.

STILLE AB cannot accept any liability for immediate or consequential damages caused by incorrect use. Use of a device for a task other than what it is intended for may damage the device and void the product warranty.

Handle the device carefully. This applies especially to the handle.

The device's performance characteristics is designed to be in line with its intended use. More detailed information regarding the device can be requested from STILLE AB when needed.

Do not store the wire saws in rotated loops, as this may affect the saw's strength and function.

The device is intended for temporary use (<60 min. under normal conditions) according to MDD 93/42/EEC, appendix 9, rule 6 (according to MDR 2017/745 appendix 8, chapter 3, rule 6).

If the handle is repaired by any service supplier not authorized by STILLE AB, all warranties are void.

CONTRAINDICATIONS

Usage outside of the intended use. Usage of the wire saw for sawing bone tissue.

CLEANING

The handle can be cleaned by soaking or applying plentifully of ethanol (70%) on device, covering all surfaces. Actuate the joint to expose all areas to the ethanol.

Leave to dry on a clean lint-free cloth or similar.

INSPECTION

Prior to each use, examine all parts for cleanliness, breaks, cracks, discoloration, corrosion, deformations and malfunction. Handles that are worn out, corroded, deformed or damaged in any other way must be sorted out.

FUNCTION TEST


Regularly examine and check the function of the handle before using it.


APPLYING

1. Pull the stockinette into place.
2. Apply felt bandage as shown.
3. Attach the self-adhesive plastic profile to the limb making sure that the profile is offset towards the bend and follows as soft a curve as possible.

Note: Do not place behind or on the ankle or other protruding bone parts. If required, attach a second profile to the other side of the limb.

4. Apply the plaster cast as usual.
5. Cut the plastic profiles to length making sure that enough remains visible for easy wire saw insertion.

REMOVAL

1. Insert the wire saw through the plastic profile.
2. Attach a handle to each end of the wire saw.
3. Start the sawing motion by carefully pulling the handles back and forth as shown, starting parallel to the limb and gradually increasing the sawing angle to a maximum of 45°. If the saw sticks, lower the angle of cutting motion back towards parallel to the limb and continue as before.


NOTE! If the patient complains of high heat generation during sawing, simply stop sawing for a moment and resume when the heat has dissipated. In case two plastic profiles are used, repeat the above for the second profile.

4. When both sides have been sawn through, remove the plaster cast by lifting the loose lid and cutting the stockinette and felt bandage. Lift the limb from the remains of the plaster cast.

PLACING WIRE SAW IN THE HANDLE

Make sure that the wire saw is mounted correctly in the handle. Correct positioning can be seen in the upper picture.

Do NOT bend the wire around one of the handle parts, as this may seriously damage the handle.

STORAGE

Devices should be stored in a clean, dry and moisture free area. Make sure that no chemicals are close to, or in the storage area.

MATERIAL

The wire saw is made of black-oxidated stainless steel W1.4528. The soft profile is made of PVC compound GFM/3E-. See separate safety data sheet. The handle is made of stainless steel W1.4021/W1.4112.

WARRANTY

STILLE AB warrants to the purchaser that its products shall be free from defects in material and workmanship. STILLE AB will repair or replace a damaged product free of charge if the damage is found to have been caused by a defect in material or workmanship. The warranty does not include normal wear and tear. The warranty is subject to the conditions that the product is used for the purpose intended, is cared for and stored according to our printed instructions and that any claims are made in writing within 30 days of the date the damage occurred. If the handles are repaired by any service supplier not authorized by STILLE AB, all warranties are void.

RETURN AND DISPOSAL OF PRODUCT

If returned for service/repair, the handle must be cleaned according to these instructions before sending back to STILLE AB or its authorized service partner. Pack the handles safely to protect them from any damage.

For disposal of defect or out-of-date product, clean it according to these instructions. The material in the handle and wire saw is recyclable as metal. The material in the plastic profile should, if possible, be sorted out as plastic.

EXPLANATION OF USED SYMBOLS

Catalogue number

Batch code

Manufacturer

Date of manufacture

Refer to instructions for use

non-sterile product

CE-mark

Do not re-use

Remaining risks weighed against the benefits of the products have been treated and minimized in the risk management handling according to EN ISO 14971. New possible risks are continuously evaluated in the PMS process.

APPLICABLE STANDARDS

EN ISO 17664

Processing of healthcare products – Information to be provided by the medical device manufacturer.

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