



Surgical perfection. For life.



STILLE

Ergonomic SuperCut Scissors

More comfort, same precision

Reducing fatigue, enhancing comfort

Ergonomic design that maintains your performance during time-consuming procedures, without compromising precision.

Surgical scissors are normally held with precision grip using the thumb and the long or ring finger. When using scissors of a traditional design, the hand usually ends up in a strenuous position which may cause tension or fatigue in the surgeon's hand after sustained and repeated use.

In order to achieve a grip that more closely resembles the hand's normal resting position, Stille has developed a range of scissors with increased space between the finger rings compared to traditional scissors. Stille's Ergo scissors have the advantage to avoid fatigue in the surgeon's hand during timeconsuming procedures. As a result, the hand joint is painless and the risk for arthritis is minimized. The dissection work becomes easier and with less effort.

STILLE-Ekberg Ergo Scissors

The portfolio of STILLE-Ekberg Ergo Scissors is developed in close collaboration with Tomas Ekberg MD, PhD, Uppsala University Hospital, Sweden. Dr. Ekberg realized the need of the new design for better balance and reduced fatigue when performing time consuming and advanced procedures.

Dr. Ekberg is a Senior Consultant and head of the department of Otorhinolaryngology and Head and Neck Surgery at the Uppsala University Hospital, Sweden.



Art. Number: 101-8019-2



Art. Number: 101-8060-2



Art. Number: 101-8175-153U-C



Art. Number: 101-8167-2-C



Art. Number: 101-8167-153U-C



Art. Number: 101-8431-172

Product Range of Ergonomic SuperCut Scissor Models

ART. NUMBER	DESCRIPTION	BLADES
101-8019-2	STILLE SuperCut Ergonomic Blepharoplasty Scissors 9.5 cm 3.75"	Curved
101-8031-2	Ragnell STILLE SuperCut Ergonomic Blepharoplasty Scissors 9.5 cm 3.75"	Curved
101-8060-2	Jabaley STILLE SuperCut Ergonomic Dissecting Scissors 12.5 cm 5"	Straight
101-8062-2	Jabaley STILLE SuperCut Ergonomic Dissecting Scissors 12.5 cm 5"	Curved
101-8163-2	Baby Metzenbaum STILLE SuperCut Ergonomic Dissecting Scissors 12.5 cm 5"	Curved
101-8173-153D	STILLE-Ekberg Metzenbaum Delicate SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved down
101-8173-153U	STILLE-Ekberg Metzenbaum Delicate SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved up
101-8173-153D-C	STILLE-Ekberg Metzenbaum Delicate Diamond SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved down
101-8173-153U-C	STILLE-Ekberg Metzenbaum Delicate Diamond SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved up
101-8175-152	Metzenbaum STILLE SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved
101-8175-182	Metzenbaum STILLE SuperCut Ergonomic Dissecting Scissors 18 cm 7"	Curved
101-8175-153D	STILLE-Ekberg Metzenbaum SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved down
101-8175-153U	STILLE-Ekberg Metzenbaum SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved up
101-8175-153D-C	STILLE-Ekberg Metzenbaum Diamond SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved down
101-8175-153U-C	STILLE-Ekberg Metzenbaum Diamond SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved up
101-8380-153D	STILLE-Ekberg Ragnell SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved down
101-8380-153U	STILLE-Ekberg Ragnell SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved up
101-8380-153D-C	STILLE-Ekberg Ragnell Diamond SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved down
101-8380-153U-C	STILLE-Ekberg Ragnell Diamond SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved up
101-8380-122	Ragnell STILLE SuperCut Ergonomic Dissecting Scissors 12.5 cm 5"	Curved
101-8380-152	Ragnell STILLE SuperCut Ergonomic Dissecting Scissors 15 cm 6"	Curved

ART. NUMBER	DESCRIPTION	BLADES
101-8166-2	Stevens-STILLE SuperCut Ergonomic Tenotomy Scissors 12.5 cm 5"	Straight
101-8167-2	Stevens-STILLE SuperCut Ergonomic Tenotomy Scissors 12.5 cm 5"	Curved
101-8167-152	Stevens-STILLE SuperCut Ergonomic Tenotomy Scissors 15 cm 6"	Curved
101-8167-182	Stevens-STILLE SuperCut Ergonomic Tenotomy Scissors 18 cm 7"	Curved
101-8167-2-C	Stevens-STILLE Diamond SuperCut Ergonomic Tenotomy Scissors 12.5 cm 5"	Curved
101-8167-153D	STILLE-Ekberg Stevens SuperCut Ergonomic Tenotomy Scissors 15 cm 6"	Curved down
101-8167-153U	STILLE-Ekberg Stevens SuperCut Ergonomic Tenotomy Scissors 15 cm 6"	Curved up
101-8167-153D-C	STILLE-Ekberg Stevens Diamond SuperCut Ergonomic Tenotomy Scissors 15 cm 6"	Curved down
101-8167-153U-C	STILLE-Ekberg Stevens Diamond SuperCut Ergonomic Tenotomy Scissors 15 cm 6"	Curved up
101-8430-172	Gorney STILLE SuperCut Ergonomic Face-lift Scissors 17 cm 6.75"	Straight
101-8431-172	Gorney STILLE SuperCut Ergonomic Face-lift Scissors 17 cm 6.75"	Curved
101-8433-202	Gorney STILLE SuperCut Ergonomic Face-lift Scissors 20 cm 8"	Curved
101-8435-202	Gorney STILLE SuperCut Ergonomic Face-lift Scissors Serrated 20 cm 8"	Curved
101-8458-172	Matarasso STILLE SuperCut Ergonomic Face-lift Scissors 17 cm 6.75"	Straight
101-8458-202	Matarasso STILLE SuperCut Ergonomic Face-lift Scissors 20 cm 8"	Curved

**Do you want to
learn more?**

Scan the QR code to
visit our website



950-G101-0003-EN-02 © 2026-01-28 STILLE. All rights reserved.

Surgical perfection. For life.



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30 Torshälla, Sweden
☎ +46 8 588 580 00
✉ info@stille.se
🌐 www.stille.se

STILLE SURGICAL INC.

840 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614

DISTRIBUTED BY: