



# STILLE imagiQ2

Upgrades & accessories



Surgical perfection. For life.

# imagiQ2 standard configuration



## Standard configuration

The STILLE imagiQ2 comes with the following configuration:

- 225 Kg / 496 lb patient weight capacity
- Black ergonomic transport handles
- Battery life of up to 2 days in normal use
- 4 swiveling wheels
- Pan handle
- Hand control + Emergency hand control
- Hand control extension cable
- Low attenuation mattress
- Fixed side rails

# imagiQ2 upgrades to fit your needs



## Mobility enhancement pack

When the table is used for frequent intra-clinic transportation, use the mobility enhancement pack for improved mobility. All wheels are designed to swivel and an additional fifth wheel is included to provide optimal mobility.

Configuration:

- 4 swiveling wheels and our proven linear 5th wheel.

Upgrade Ref: 535-1711

# imagiQ2 upgrades to fit your needs



## Load upgrade pack

Maximize the capabilities of the imagiQ2 with increased patient weight capacity up to 661 lb / 300 kg. This premium upgrade has silver ergonomical transport handles.

Please note that this option must be specified on original order before delivery. It can't be retrofitted to an already delivered unit.

Upgrade Ref: 535-1710

# imagiQ2 upgrades to fit your needs



## Performance pack

For maximum uptime, choose the performance pack. It includes a high-performance long-term battery pack with up to 1 week of durability.

This upgrade can be chosen at any time.

Upgrade Ref: 535-1712

# imagiQ2 accessory highlights



**Extra pan handle**  
Extra pan handle allows panning the table from different sides



**Detachable short side rail**  
Short 8 in / 20 cm rail system enabling accessory mounting and 360° imaging capabilities.



**Translucent fistula arm board**  
Provides an unobstructed image and is connected directly to the carbon fiber top.



**Adjustable procedure headrest**  
Fully translucent headrest, ensures easy intubation and superior carotid visualization. It fits the free end of the table and is side rail mounted.



**imagiQ2 radiation shield**  
Tailored and easy to fit protection against scattered radiation with lead rubber (0.5 Pb) material.



**Detachable side rail**  
One third rail system enabling full table length side rails, accessory mounting and optional 360° imaging capabilities.



**Anesthesia screen**  
Provides a sterile barrier separating the surgical field from the patient's head.



**Catheter tray**  
Provides an extended work space, which is attached directly to the table. Available in three lengths.



**Adjustable arm board**  
Nontranslucent board provides ideal support for the patient's arm during IV.



**imagiQ2 foot control**  
Compact multipedal footswitch for discrete placement beneath the table with soft touch buttons and a light operating force.

The imagiQ2 offers a large line of accessories for vascular tables and can be upgraded with features and detachable table tops for different purposes when it is needed.

# imagiQ2 accessory catalogue


Accessory	Part number	Standard	Option	Market availability
Pan handle DIN	535-1760	x		All excl US
Pan handle US	535-1761	x		US only
Detachable side rail US	535-1720		x	All
Detachable side rail DIN	535-1721		x	All
Detachable side rail UK	535-1722		x	All
Short detachable 8 inch rail US*	535-1724		x	All
Short detachable 8 inch rail DIN*	535-1725		x	All
Short detachable 8 inch rail UK*	535-1726		x	All
Fistula arm board	535-1730		x	All
Head rest, adjustable	535-1740		x	All
Head rest fixed	535-1741		x	All
Catheter tray 1000 mm	535-1750		x	All
Catheter tray 1500 mm	535-1751		x	All
Catheter tray 500 mm	535-1752		x	All
Foot Control	535-1762		x	All
Adjustable armboard DIN	514-50-3SR025		x	All excl US
IV pole	514-100401		x	All excl US
Anesthesia frame	514-100311		x	All excl US
Basic armboard	512-5000		x	US only
Anesthesia screen	512-5030		x	US only
IV pole w/o clamp	512-5010		x	US only
IV pole clamp	512-5013		x	US only
Footboard carbon fiber with pad	512-5040		x	US only
Deluxe maximum restraint	O-RSX2		x	US only
Radiation shield base kit grey**	535-1770		x	EU only
Radiation shield base kit blue**	535-1771		x	EU only
Radiation shield loop grey. Extra shield**	535-1772		x	EU only
Radiation shield loop blue. Extra shield**	535-1773		x	EU only
Radiation shield base kit grey**	535-1770-US		x	US only
Radiation shield base kit blue**	535-1771-US		x	US only
Radiation shield loop grey. Extra shield**	535-1772-US		x	US only
Radiation shield loop blue. Extra shield**	535-1773-US		x	US only

\* The new detachable 20 cm / 8 in side rail is not to be used with a retractor system, catheter trays 1000 mm and 1500 mm

\*\*The system is designed to fit the imagiQ2 and consists of 1 mounting piece that is placed under the mattress and 3 radiation shield curtains. Extra shield curtains, which also can be fitted to the side rail if customer only requires a single curtain and not full whole system

## imagiQ2 upgrades

Upgrade features	Description
Mobility enhancement pack	Ref: 535-1711
Load upgrade pack	Ref: 535-1710
Performance pack	Ref: 535-1712



## STILLE - Surgical perfection. For life.

STILLE is one of the oldest medical device companies in the world. At STILLE, we develop, produce and sell high-quality surgical instruments and surgical imaging tables to physicians worldwide, and has done so for 170 years.

STILLE's offering consists of a strong brand and products renowned for their superior quality. We manage our business in two distinct business areas—Surgical Instruments and Surgical imaging tables.

### Surgical instruments

In our Surgical Instruments business area, we have a complete offering of handcrafted, forged instrumentation for all fields of open surgery, with a special emphasis on cardiovascular surgery, and reconstructive and aesthetic plastic surgery.

Our instruments feature delicate design, industry-leading durability and a unique feel, which has made STILLE a well-known and highly respected supplier by surgeons throughout the world.

### Surgical imaging tables

In our Surgical imaging tables business area, we offer specialty tables for minimally invasive procedures, with pricing, safety and features optimized for both in- and outpatient settings. The STILLE imagiQ2 surgical imaging table with its low-dose capabilities, is optimized for advanced mobile Hybrid OR settings as well as peripheral and interventional procedures.

Learn more at [www.stille.se](http://www.stille.se)



Surgical perfection. For life.

#### STILLE AB

Gustav III:s Boulevard 42  
SE-169 73 Solna, Sweden  
Phone: +46 8 588 580 00  
Fax: +46 8 588 580 05  
[info@stille.se](mailto:info@stille.se)

[www.stille.se](http://www.stille.se)

#### STILLE Surgical Inc.

943 Parkview Blvd  
Lombard, IL 60148  
United States  
Phone: 1.800.6651614

#### Distributed by: