


Low-profile applicators and removers to suit every surgical need.

Like our T2 Clips, Mizuho applicators and removers feature a slim, low-profile design for optimal visibility during placement or removal.



Diverse Instrument Selection

Mizuho offers a full variety of applicators to complement T2 Clips. There are Short (70 mm), Standard (90 mm), and Long (110 mm) lengths available, as well as straight and angled shapes. The slim, ultra-light, perfectly balanced applicators feature meticulously milled jaw grooves to facilitate a sure, no-slip grip and vibration-free release.

Uniquely Color Coded for Safety, Ease of Use



T2 Clips are fully compatible with Mizuho's familiar line of Standard and Mini Applicators. By color coordinating the applicator's jaw with the color of the clip spring, the correct applicator/clip combination can be verified at a glance, ensuring proper instrument selection during surgery.



SANO Multi-Angle Applicator



SANO Double-Action Applicator



Ordering information.

Aneurysm Clip Cases

Product Code	Product Name
17-010-50	Sugita Titanium 2 Clip Sterilization Case (Full set)
17-010-50S	Sugita Titanium 2 Clip Sterilization Case (Standard/Large clips)
17-010-50T	Sugita Titanium 2 Clip Sterilization Case (Mini/Temporary/Standard clips)
17-010-50F	Sugita Titanium 2 Clip Sterilization Case (Fenestrated/Large clips)
17-010-50D	Sugita Titanium 2 Clip Sterilization Case (Defected Fenestrated/Long clips)

Standard Clips

Product Code	Description
17-001-01	T2 8.5 Straight 7mm
17-001-02	T2 10.5 Straight 10mm
17-001-03	T2 11.0 Bayonet 7mm
17-001-04	T2 9.5 Bent 8mm
17-001-05	T2 11.0 Sideward Bent 8mm
17-001-06	T2 10.0 Sideward Curved 9mm
17-001-07	T2 12.0 Sideward Curved Bayonet 8.5mm
17-001-08	T2 9.5 Curved 9mm
17-001-09	T2 10.0 Slightly Curved 9mm
17-001-10	T2 11.5 Straight 12mm
17-001-11	T2 13.5 Sideward Bent 12mm
17-001-12	T2 11.0 Bent 12mm
17-001-13	T2 13.0 Bayonet 10mm
17-001-14	T2 14.0 Bayonet 12mm
17-001-15	T2 13.5 Straight 15mm
17-001-16	T2 14.0 Slightly Curved Bayonet 12mm
17-001-17	T2 15.5 Slightly Curved 18mm
17-001-18	T2 15.5 Straight 18mm
17-001-20	T2 8.0 L-Shaped 5mm
17-001-21	T2 8.0 L-Shaped 7.5mm
17-001-22	T2 8.0 L-Shaped 10.0mm
17-001-23	T2 9.0 J-Shaped 6mm
17-001-24	T2 9.0 J-Shaped 7.5 mm
17-001-25	T2 9.5 J-Shaped 9.0mm
17-001-47	T2 10.5 45° Angle 5mm
17-001-48	T2 11.5 45° Angle 7.5mm
17-001-49	T2 12.5 45° Angle 10mm
17-001-50	T2 13.0 Curved 14mm
17-001-60	T2 8.0 L-Shape w/ Curve 7.5mm
17-001-61	T2 8.0 L-Shaped w/ Curve 10mm
17-001-62	T2 9.0 1/4 Curved 7.5mm
17-001-63	T2 9.5 1/4 Curved 8.6mm
17-001-64	T2 10.5 Slightly Bent 8mm
17-001-65	T2 11.0 Slightly Curved 11mm
17-001-66	T2 11.0 Curved 11mm
17-001-68	T2 10.0 Slightly Curved 9mm
17-001-69	T2 9.5 Curved 9mm
17-001-70	T2 17.0 Straight 21mm (19A)
17-001-73	T2 11.0 1/4 Curved 11mm

Standard Temporary Clips

Product Code	Description
17-001-51	T2 Temp. 9.0 Straight 7mm
17-001-52	T2 Temp. 10.5 Straight 10mm
17-001-53	T2 Temp. 11.5 Bayonet 7mm
17-001-54	T2 Temp. 9.5 Bent 8mm
17-001-55	T2 Temp. 11.0 Slightly Curved 11mm

Standard Fenestrated Clips

Product Code	Description
17-001-26	T2 11.0 Fen. 3.5 Straight 6mm
17-001-27	T2 13.0 Fen. 3.5 Straight 9mm
17-001-28	T2 14.5 Fen. 3.5 Straight 12mm
17-001-29	T2 12.0 Fen. 5.0 Straight 6mm
17-001-30	T2 14.0 Fen. 5.0 Straight 9mm
17-001-31	T2 15.5 Fen. 5.0 Straight 12mm
17-001-32	T2 8.5 Fen. 3.5 Bent 5mm
17-001-33	T2 9.0 Fen. 3.5 Bent 7.5mm
17-001-34	T2 10.0 Fen. 3.5 Bent 10mm
17-001-35	T2 9.5 Fen. 5.0 Bent 5mm
17-001-36	T2 10.0 Fen. 5.0 Bent 7.5mm
17-001-37	T2 11.0 Fen. 5.0 Bent 10mm
17-001-38	T2 7.5 Fen. 3.5 L-Shaped 5mm
17-001-39	T2 7.5 Fen. 3.5 L-Shaped 7.5mm
17-001-40	T2 7.5 Fen. 3.5 L-Shaped 10mm
17-001-41	T2 8.5 Fen. 5.0 L-Shaped 5mm
17-001-42	T2 8.5 Fen. 5.0 L-Shaped 7.5mm
17-001-43	T2 8.5 Fen. 5.0 L-Shaped 10mm
17-001-44	T2 10.5 Fen. 5.0 Bayonet 5mm
17-001-45	T2 12.0 Fen. 5.0 Bayonet 7.5mm
17-001-46	T2 13.5 Fen. 5.0 Bayonet 10mm
17-001-74	T2 8.5 Fen. 5.0 L-Shaped Cvd. Left 5mm
17-001-75	T2 8.5 Fen. 5.0 L-Shaped Cvd. Left 7.5mm
17-001-76	T2 8.5 Fen. 5.0 L-Shaped Cvd. Left 10mm
17-001-77	T2 8.5 Fen. 5.0 L-Shaped Cvd. Right 5mm
17-001-78	T2 8.5 Fen. 5.0 L-Shaped Cvd. Right 7.5mm
17-001-79	T2 8.5 Fen. 5.0 L-Shaped Cvd. Right 10mm
17-002-04	T2 11.0 Fujita Fen. 5.0 Angled 5mm
17-002-05	T2 12.0 Fujita Fen. 5.0 Angled 7.5mm
17-002-07	T2 9.0 Fujita Fen. 5.0 70° Angled 5mm (R)
17-002-08	T2 10.0 Fujita Fen. 5.0 70° Angled 7.5mm (R)
17-002-09	T2 9.0 Fujita Fen. 5.0 70° Angled 5mm (L)
17-002-10	T2 10.0 Fujita Fen. 5.0 70° Angled 7.5mm (L)
17-002-11	T2 8.0 Fujita Fen. 5.0 L-Shaped Angled 5mm (R)
17-002-12	T2 8.0 Fujita Fen. 5.0 L-Shaped Angled 7.5mm (R)
17-002-13	T2 8.0 Fujita Fen. 5.0 L-Shaped Angled 5mm (L)
17-002-14	T2 8.0 Fujita Fen. 5.0 L-Shaped Angled 7.5mm (L)

Mini Clips

Product Code	Description
17-001-80	T2 Mini 5.5 Straight 4mm
17-001-81	T2 Mini 7.0 Straight 6mm
17-001-83	T2 Mini 6.5 Curved 5mm
17-001-84	T2 Mini 7.0 Bent 5mm
17-001-85	T2 Mini 5.5 Curved 4mm
17-001-86	T2 Mini 7.5 Sideward Curved 7mm
17-001-87	T2 Mini 6.0 Curved 5mm
17-001-88	T2 Mini 7.5 Curved 7mm
17-001-89	T2 Mini 8.5 Bayonet 5mm
17-001-95	T2 Mini 5.5 Straight Slim Tip 4mm
17-001-96	T2 Mini 7.0 Straight Slim Tip 6mm
17-001-97	T2 Mini 7.0 Curved Slim Tip 6.7mm
17-001-98	T2 Mini 8.0 Curved Slim Tip 8.2mm
17-001-99	T2 Mini 6.5 Curved Slim Tip 5.2mm
17-002-15	T2 Mini 8.0 Straight 8mm
17-002-16	T2 Mini 7.0 Curved 6.5mm
17-002-17	T2 Mini 6.5 Slightly Curved 5.5mm
17-002-18	T2 Mini 7.0 Slightly Curved 6.5mm
17-002-19	T2 Mini 8.5 Sideward Bent 6mm
17-002-20	T2 Mini 10.0 Bayonet 7mm
17-002-21	T2 Mini 6.5 Curved Slim Tip 5.5mm

Mini Temporary Clips

Product Code	Description
17-001-56	T2 Mini Temp. 7.0 Straight 6mm
17-001-57	T2 Mini Temp. 7.0 Bent 5mm
17-001-58	T2 Mini Temp. 7.0 Straight Slim Tip 6mm

Sugita Titanium 2 Aneurysm Clip Applicators and Removers

Product Code	Description
17-012-52	Standard Straight 90mm
17-012-54	Standard 15° Angled 90mm
17-012-55	Standard 30° Angled 90mm
17-012-56	Standard Remover Straight 90mm
17-013-52	Mini Straight 90mm
17-013-54	Mini 15° Angled 90mm
17-013-55	Mini 30° Angled 90mm
17-013-56	Mini Remover Straight 90mm
17-012-51	Standard Straight 70mm
17-012-53	Standard Straight 110mm
17-012-57	Standard Remover 30° Angled 90mm
17-012-58	Standard Laterally Angled 110mm
17-012-09	T2 Sano Double Action Applier for Std. Clips
17-012-10	T2 Sano Multi Angle Applier for Std. Clips
17-013-51	Mini Straight 70mm
17-013-53	Mini Straight 110mm
17-013-57	Mini Remover 30° Angled 90mm
17-013-58	Mini Laterally Angled 110mm
17-013-09	T2 Sano Double Action Applier for Mini
17-013-10	T2 Sano Multi Angle applier for Mini