



Stille Forceps for  
Plastic Surgery

## Stille - Surgical Perfection

Stille is one of the oldest medical device companies in the world. We develop, manufacture and sell proprietary high-quality surgical instruments and specialty tables to leading surgeons worldwide. We have done so for more than 170 years.

Stille's offering consists of a strong brand and products renowned for their superior quality. We manage our business in two distinct business areas - Surgery and Patient Positioning.

Stille surgical instruments are considered state of the art by many surgeons around the world. All Stille instruments are handcrafted as an art by highly skilled instrument makers in Eskilstuna, Sweden.

The slender design and excellent balance are the characteristics which make the Stille Instruments easily recognizable. Their highly polished surface and delicate design make Stille instruments lighter and easier to keep clean, thus more durable against discoloration and corrosion. All this provides Stille instruments with a unique and outstanding feeling.



Surgical Perfection. For life.

**Stille AB**  
Gustav III:s Boulevard 42  
SE-169 73 Solna, Sweden  
Phone: +46 8 588 580 00  
Fax: +46 8 588 580 05  
sales@stille.se

**Stille Surgical Inc.**  
943 Parkview Blvd  
Lombard, IL 60148  
United States  
Phone: 1.800.6651614  
sales.us@stille.se

Distributed by:

[www.stille.se](http://www.stille.se)

Art. no. 950-102-EN13-01 © 2013-02-06 Stille. All rights reserved.



## Stille Forceps for Plastic Surgery


Surgical precision all the way to the tips



Surgical Perfection. For life.



# Perfect tissue handling. Every time.



102-9077-20  
**Crafoord-Stille Dressing Forceps**  
Serrated, 20 cm 8"

As an extension of the surgeon's thumb and forefinger, forceps need to grasp and release with the same sensitivity as a pair of fingertips.

102-9025-1  
**McIndoe Stille Dressing Forceps**  
Diamond Jaws  
15.5 cm 6.25"

102-9022  
**Adson Stille Tissue Forceps**  
Cross-serrated, 1x2 teeth  
12 cm 4.75"

102-9252-25  
**Waugh-Brophy Stille Tissue Forceps**  
1x2 teeth, Cross-serrated, 25 cm 10"

As an extension of the plastic surgeon's thumb and forefinger, forceps need to grasp and release with the same sensitivity as a pair of fingertips. The claws on surgical forceps are specially designed to make a firm, yet gentle, grasp of the tissue.

The anatomical forceps are designed, as much as possible, to take an atraumatic grip on sensitive tissue. It is the fine tuning that creates the perfect balance with a smooth, resilient spring action and a firm grip. The surface grooves increase gripping power even further. This gives the plastic surgeon full control while causing as little harm to the tissue as possible.

Forceps are among the surgical instruments that founder Albert Stille learned to make from scratch during his apprenticeship in the 1830's. Forceps are seemingly simple instrument but still demands ultra-precision in the manufacturing process. Over the years, Stille's forceps continue to be refined and are today recognized for their unique feel and excellent balance in combination with unsurpassed longevity.

102-9023-12  
**Adson-Brown Stille Tissue Forceps**  
Side-grasping teeth  
12 cm 4.75"

102-9027  
**Gillies Stille Tissue Forceps**  
Serrated, 1x2 fine teeth  
15.5 cm 6.25"

102-9021-1  
**Adson-Callison Stille Tissue Forceps**  
Delicate tying platform  
1x2 teeth, 12 cm 4.75"

102-9020  
**Adson Stille Dressing Forceps**  
Serrated, 12 cm 4.75"

102-9852-20  
**Waugh-Brophy Stille Diathermy Forceps**  
Insulated, 1x2 teeth, Cross-serrated, 20 cm 8"

102-9845-20  
**Björk-Stille Diathermy Forceps**  
Insulated, 20 cm 8"