




STILLE Surgical Instruments

Quality tools behind the talent





"The world's best instruments,
for the world's best surgeons"

STILLE Surgical Instruments

For over 175 years, we have developed and manufactured the best surgical instruments for the world's most demanding surgeons. We would like to extend a sincere thank you to all our current customers and welcome our future customers to experience the STILLE difference. In this catalog we present a selection of STILLE original surgical instruments.

Precision, durability and feel are characteristic qualities of all STILLE instruments. The majority are handcrafted by our highly skilled artisans in Eskilstuna, Sweden. Our instrument manufacturing process from melting steel rods to finished product is a long one, involving multiple stages of development. Our method of crafting STILLE instruments involves only the highest-grade steel, which provides a unique feel and reliability to the surgeon in any surgical specialty.

In this catalogue you will find a wide selection of instruments, including the original supercut scissor developed in 1984 and our new Diamond supercut design. Quality goes into our STILLE forceps, needle holders, clamps, rongeurs, curettes and retractors, as well as, a fine selection of stainless and titanium micro instruments. Common to all of our instruments is that they are smooth and consistently perform to your expectations. Our outstanding balance, solid design and ergonomics provides the surgeon with a reliable instrument throughout long and complex procedures.

We at STILLE are proud to have developed the world's first double-action bone ronguer. We also developed the world's first osteotome with the classic 8-sided handle. Many STILLE instruments have been copied, but a copy cannot measure up to the original.



Dr. Mark Jewell

"My patients deserve nothing but the best.
STILLE has been the natural extension of my
hands for over 30 years"

Mark Jewell, M.D., P.C.
Plastic Surgeon
www.markjewellmd.com

Scissors



STILLE Eye Scissors

Delicate eye scissors with sharp or blunt tips and large finger rings. These are also called iris scissors and are used for nerves, tendons and vessels.

101-8003



STILLE SuperCut XE Eye Scissors

SuperCut XE eye scissors with sharp or blunt tips and extremely sharp edges. Beveled blades.

101-8015-1



Kaye STILLE SuperCut XE Dissecting Scissors

SuperCut XE version for improved cutting performance.

101-8021-1



Jabaley STILLE SuperCut XE Ergonomic Dissecting Scissors

Ergonomic model of Super Cut XE scissors for dissection of tough tissue in hand surgery. Beveled blades.

101-8060-2



**Baby Metzenbaum
STILLE SuperCut XE
Dissecting Scissors**

Super Cut dissecting scissors for tough tissue. Used in delicate surgery requiring small instruments.
101-8161



**Mayo STILLE
Dissecting Scissors**

The most used scissors, for all types of surgery. Used for routine dissection and cutting in tissue. Beveled blades.
101-8080-17



**Metzenbaum STILLE
SuperCut XE Dissecting
Scissors**

The popular Metzenbaum dissecting scissors with Super Cut edges. Rounded/beveled blades.
101-8175-18



**Metzenbaum STILLE
Diamond SuperCut**

Popular dissecting scissors with Super Cut edges. Rounded/beveled blades. Extremely hard and low-reflecting diamond coating.
101-8175-18-C



**Ragnell STILLE
Dissecting Scissors**

Dissecting scissors with flat, rounded tips. Often used in plastic surgery. Is also called Killner by some manufacturers.
101-8380-15

Scissors



Stevens-STILLE SuperCut XE Tenotomy Scissors

Very popular SuperCut XE scissors for precise cutting and dissection of fine tissues. Rounded blades.
101-8167-151



Peck-Joseph STILLE Nasal Scissors

Double-beveled scissors for rhinoplasty.
101-8125



STILLE SuperCut XE Vascular Scissors

Super Cut scissors with angled blades for better reach of vessels etc. Used in vascular surgery.
101-8200-60



Matarasso STILLE SuperCut XE Face-lift Scissors

SuperCut XE face-lift scissors with double-beveled semi-sharp outer blades.
101-8458-20



STILLE Suture Scissors

Suture removal scissors.
101-8584

Micro Instruments



STILLE SuperCut Micro Scissors

Fine micro scissors with straight, sharp tips or curved blunt tips. All models with SuperCut blades.
117-A110-18



STILLE SuperCut Micro Vascular Scissors

SuperCut XE micro scissors available in five different blade angles for fine vascular surgery.
117-A120-45



STILLE Micro Forceps

Well balanced and ergonomic micro forceps of ultra high quality. With tying platform and straight or curved jaws.
117-B110-18



STILLE Micro Ring Forceps

Round handle, with fine ring tips, counter-balanced.
117-B140-18

Micro Instruments



STILLE Micro Tissue Forceps

Round handle, with
1x2 fine teeth.

117-B150-15



STILLE Micro Needle Holder

Same fine micro needle
holder, with lock.

117-C121-18



STILLE Titanium Micro Forceps

Round handle, with tying platform,
straight or curved tips.

117-B120-18-TI



STILLE Titanium Micro Forceps

Round handle, with tying platform,
counter-balanced, straight or
curved tips.

117-B130-21-TI



STILLE Titanium Micro Ring Forceps

Round handle, with fine ring tips,
counter-balanced.

117-B140-21-TI



STILLE Titanium Micro DeBakey Forceps

Round handle, with 1.2 mm
fine atraumatic DeBakey tips.

117-B160-18-TI



STILLE Titanium Micro Needle Holder

Round handle, fine 0.8 mm Diamond
Jaws and Regular Box Lock for sutures
7-0 and smaller. With or without ratchet.

117-C110-18-TI



STILLE Titanium Micro Needle Holder

Round handle, 1.2 mm Diamond Jaws
and Streamline Box Lock for sutures 5-0
and smaller. With or without ratchet.

117-C125-21-TI

Forceps



Adson STILLE Tissue Forceps

Tissue forceps with very delicate points and 1x2 fine teeth. Used in most specialties for delicate work.

102-9022



Gillies STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth and stabilizing pin. An excellent forceps for the surgeon who wants an Adson-size tip, but with greater length.

102-9027



STILLE De Bakey Vascular Forceps

Atraumatic forceps with De Bakey jaws for handling of blood vessels and other delicate tissue. Can be insulated on request.

102-9312



Waugh-Brophy STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth. Very fine tips, like Adson. Useful when long and slender instruments are needed.

102-9252-25



McIndoe STILLE Dressing Forceps

Dressing forceps with stabilizing pin and delicate jaws for grasping fine tissue.


102-9025



STILLE Dressing Forceps

Dressing forceps with square shaped and serrated points used in all surgical specialties.

102-9030-15



Björk-STILLE Diathermy Forceps

Insulated forceps for monopolar diathermy. Serrations allow grasping of delicate structures.


102-9845-20



Crafoord-STILLE Dressing Forceps

Dressing forceps with rounded points. The perfect balance and delicate grip gives this instrument a unique feeling.

102-9077-20



Olivecrona-STILLE Dura Forceps

Dura forceps with 2x3 fine teeth. The nonpenetrating tips give the surgeon maximum grip without damaging the tissue.

102-9244-15



STILLE Tissue Forceps

Tissue forceps with 2x3 teeth. Used when more grasping power than 1x2 teeth is needed.

102-9062-20

Needle Holders



**Derf STILLE
Needle Holder**

Needle holder with rather short TC jaws, popular in plastic surgery. Intended for 3/0 suture and finer.

104-15542



**Halsey STILLE
Needle Holder**

Needle holder with TC jaws used in many procedures. One model with ergonomically designed finger ring for thumb.

104-15560



**Olsen-Hegar STILLE
Needle Holder**

Popular needle holder with scissors for suturing and cutting with the same instrument. Intended for 0 suture and finer.

104-15590-131



**Webster/Kleinert STILLE
Needle Holder**

Delicate needle holder for fine needles. With smooth jaws or TC jaws.


104-15559



**Mayo-Hegar STILLE
Needle Holder**

Commonly used general purpose needle holder with TC jaws. Intended for 0 suture and finer.

104-15561-181




STILLE Crile-Wood Needle Holder

General purpose needle holder with TC jaws. Intended for 2/0 suture and finer.
104-15565-201



Gillies STILLE Needle Holder

Needle holder with scissors. Ergonomic design on shanks and finger rings. One rather large ring handle is set at an angle for fitting the top of the thumb. Without ratchets. With or without TC jaws. The scissor part enables suturing and cutting with one instrument. Intended for 2/0 suture and finer.
104-15551



STILLE-Frimand DUO-GRIP Needle Holder


This new and innovative needle holder reduces the suturing time and tissue manipulation as well as offering better ergonomics for the surgeon.
104-15665-151



Rapp-STILLE Knot-holding Forceps

A 13 cm knot-holding forceps with smooth jaws. Ideal instrument to use as a tool for holding the suture knot in narrow areas.
104-15140

Vascular Clamps



STILLE Multi-purpose Clamp

Multi purpose clamps with angled, atraumatic DeBailey jaws.
107-10825




Satinsky STILLE Tangential Occlusion Clamp

Tangential occlusion clamps with atraumatic DeBailey jaws.
107-10841



STILLE Anastomosis Clamp

Anastomosis clamps with atraumatic DeBailey jaws.
107-10860



STILLE DeBailey Aortic-Aneurysm Clamp

Aortic aneurysm clamp with atraumatic DeBailey jaws.
107-10849



Ericsson-STILLE Carotid Clamp


Carotid clamp with angled atraumatic DeBailey jaws.
107-10880

Other Clamps



Backhaus STILLE Towel Clamps

Towel clamps in two different versions. Sharp/sharp or ball/socket for paper towels.
107-10910



Babcock STILLE Intestinal Grasping Forceps

Intestinal grasping forceps in two different lengths.
107-12370-16



Allis STILLE Intestinal Grasping Forceps


Tissue-grasping forceps for secure grip of tissue.
107-13120-15



Foerster STILLE Sponge Holding Forceps

Sponge holding forceps with several different loops and lengths.
107-12485


Bone Cutters and Rongeurs



STILLE Extraction Plier

Original STILLE double-action joint facilitates removal of pins and wires.

105-16020



STILLE Bone Cutting Forceps

Original STILLE double-action angled bone cutting forceps.

105-16811



Beyer STILLE Bone Rongeur

Bone rongeur with original STILLE double-action joint and two different jaw widths.

105-17210



Leksell STILLE Bone Rongeur

Bone rongeur with original STILLE double-action joint and sideways slightly curved jaws.


105-17260



Semb STILLE Bone Holding Forceps

Bone holding forceps with or without locking mechanism.

105-16150



STILLE Wire Cutter

Double joint wire cutter with TC-inserts. For hard wire up to 1.7 mm and soft wire up to 2.0 mm.

105-16223



Kerrison STILLE Laminectomy Punch

The popular Kerrison punches with a unique rubber coated handle for improved ergonomics. No spring between the handles improves the accessibility in cleaning and decreases the risk for getting caught with the gloves.

105-17642-5



Caspar STILLE Rongeur

Rongeur for intervertebral discs, with semi-detachable shaft. Available with or without coating. Fenestrated and serrated jaws in three different sizes.

105-13956-18

Artery Forceps



Halsted-Mosquito STILLE Artery Forceps

Hemostat without teeth, used in small operative fields. For plastic surgery, neuro, hand or any delicate surgery.
103-10032-121



Pean STILLE Artery Forceps

General utility forceps used in many surgical disciplines.
103-10015-131



Crile STILLE Artery Forceps

General purpose hemostat without teeth. Used in relatively small operating fields, when more strength is needed than a Mosquito can offer.
103-10182-161



Leriche STILLE Artery Forceps

Fairly tapered hemostat. For use in small operative fields.
103-10192

Skin Hooks & Retractors



Gillies STILLE Skin Hook

Skin hook with small or large hook on long flat shaft.

108-24110



Joseph STILLE Skin Hook

Skin hook with sharp prongs. Available in five different models. Popular in plastic surgery.


108-25012-10



STILLE Breast Retractor

Useful double-sided retractor for breast augmentation. Serrated blades in two different sizes and shapes. Ergonomic handle for maximum comfort.

108-26038



Weitlaner STILLE Retractor

General purpose self-retaining retractor used in many procedures.

108-26254

Bone Instruments



STILLE Bone Curette

Oval and sharp bone curettes with widths from 2 mm to 12 mm.

111-20502-4



Wiberg STILLE Raspatory

A straight raspatory in four different sizes.

111-22171-7



Howarth STILLE Nasal Raspatory

Double-ended nasal raspatory with sharp hooks.

111-22370



Freer STILLE Elevator

Double-ended elevator sharp/blunt with rounded, delicate handle.

111-22140-6



Fomon STILLE Rasp

Double-ended rasp with TC-inserts for better durability. Low profile for better access in narrow areas.

111-22860



STILLE Hand Osteotome

Delicate osteotome with phenolic handle. Intended for hand surgery, but also useful in other fine orthopaedic work.

111-21456-4



Lambotte STILLE Osteotome

Straight osteotome available in 6 different sizes.

111-21441-15



Lexer STILLE Osteotome

Osteotome with phenolic handle.


111-21460-15



STILLE Gouge

Original STILLE gouge. Available in five different sizes from 6 mm to 15 mm.


111-21315-8



STILLE Mallet

Stainless steel mallet.

111-21820





For a complete presentation of all STILLE instruments, see our product catalogue at; www.stille.se



CRAFTSMANSHIP ALL THE WAY. Creating quality instruments is an art form – from instrument maker to finished product. STILLE has mastered this process every step of the way. Passion, professionalism and dedication are common attributes. Through our close collaboration with surgeons all over the world, our instruments are fine-tuned and developed with the precision they consistently expect from STILLE. Our product range is specialized in all surgical specialties including plastics, cardiac and orthopedic. True craftsmanship, outstanding performance and superior quality are the unique characteristics that make STILLE the hands down choice of surgeons worldwide.

STILLE AB

Ekbacken 11
SE-644 30 Torshälla, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
info@stille.se
sales@stille.se

STILLE Surgical Inc.

915 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1-800-665-1614

www.stille.se