



Surgical perfection. For life.

STILLE Medstone CV Cardiovascular/ Interventional Radiology

High Quality Hospital Grade Surgical Imaging Tables with Free Floating or Motorized Capabilities

The STILLE Medstone c-arm tables contribute to enhancing efficiency while representing an affordable solution with quality and reliability provided by Swedish technology.



Photo shows Medstone3 CV
(500-1025)

STILLE Medstone3 CV Floating table Technical Specifications for Standard Configuration

- 36" (91,5cm) longitudinal floating travel
- 8" (20cm) lateral floating travel
- Carbon fiber table top 24"x82" (61x208cm)
- AP imaging length 65.7" (167cm);
3D imaging length 30.7" (78cm)
- Motorized height adjustment 26.6"-42.4" (67.6-107.6cm)
- Weight capacity 500lb (227kg)
- 230Vac, 2A, 50/60Hz and 120Vac, 3,6A, 50/60Hz
with 230Vac or 120Vac style line cord
- Four hospital grade locking casters
- Float control unit with magnetic locking release
(side rail mounted)
- Table mattress 2" (5cm) and 3 patient straps
- Foot control and hand control

Shortening Procedure Times and Enhancing OR Efficiency

STILLE has the longest published longitudinal float for mobile surgical imaging tables. The Medstone3 CV table with free floating capabilities allows for head to toe imaging with accuracy, speed and freedom of movements to improve efficiency in the OR.

STILLE Medstone3 CV Surgical Imaging Table

- Free floating top
(industry leading longitudinal float range)
- Motorized elevation
- Patient weight capacity 500lb (227kg)
- Carbon fiber table top with 1mm Aluminium Equivalence
- Four hospital grade locking casters
- Float control unit with magnetic locking release

Model	Medstone3 CV	Medstone5 CV
Part number	500-1025	500-1002
Table top	24"x82" (61x208cm)	22"x82" (56x208cm)
Height adjustment	26.6"-42.4" (67.6-107.6cm)	26.78"-39.76" (65.5-101.0cm)
Trendelenburg	x	± 15°
Lateral tilt	x	± 15°
Longitudinal travel	36" (91.5cm) free floating	15.75" (40cm) motorized
Lateral travel	8" (20cm) free floating	6" (15cm) motorized
Infectious control covers	x	✓
Fixed side rails (on both long table sides)	36.51" (93cm)	11,97" (30.4cm)
Accessory Package	B	A

Specification tolerance: ±1°, ±.5" (1.27cm)

A= Table pad, wired hand control

B=Same as A plus foot control, float control unit and 3 patient straps

Optional Accessories

Accessories for Medstone3 CV

Part number	Description
505-1118	Carbon fiber fistula arm board with pad for imaging and injection
505-1119	Accessory table clamp with 5.9" (15cm) rail to attach float control unit to head end of table
505-1131	Catheter tray, to be attached to side rail
505-1121	2 accessory clamps with rail to attach catheter tray to head end of table
505-1179	2 detachable side rails length 31.69" (80.5cm)
505-1112	Patient arm board with pad and clamp for attachment directly to the table top
514-100401	IV-pole with clamp (side rail mounted)
514-100311	Anesthesia frame with clamps (side rail mounted)
535-1770	X-ray protective kit grey, 6 single panels; attachment to table top with self-adhesive Velcro © strips
535-1771	X-ray protective kit blue, 6 single panels; attachment to table top with self-adhesive Velcro © strips
535-1772	X-ray protective loop grey, 5 connected panels, attachment to side rails or table top with self-adhesive Velcro © strips
535-1773	X-ray protective loop blue, 5 connected panels, attachment to side rails or table top with self-adhesive Velcro © strips

Accessories for Medstone5 CV

Part number	Description
505-1175	2 detachable side rails length 44.29" (112.5cm)
505-1112	Patient arm board with pad and clamp for attachment directly to the table top
514-100401	IV-pole with clamp (side rail mounted)
514-100311	Anesthesia frame with clamps (side rail mounted)
535-1770	X-ray protective kit grey, 6 single panels; attachment to table top with self-adhesive Velcro © strips
535-1771	X-ray protective kit blue, 6 single panels; attachment to table top with self-adhesive Velcro © strips
535-1772	X-ray protective loop grey, 5 connected panels, attachment to side rails or table top with self-adhesive Velcro © strips
535-1773	X-ray protective loop blue, 5 connected panels, attachment to side rails or table top with self-adhesive Velcro © strips

STILLE – Surgical perfection since 1841

STILLE is one of the oldest medical device companies in the world founded at Karolinska University Hospital in Sweden back in 1841. At STILLE, we develop, produce and sell high-quality surgical instruments and surgical imaging tables to physicians worldwide.

STILLE's offering consists of a strong brand and products renowned for their superior quality. We manage our business in two distinct business areas—Surgical Instruments and Surgical imaging tables.

Surgical instruments

In our Surgical Instruments business area, we have a complete offering of handcrafted, forged instrumentation for all fields of open surgery, with a special emphasis on cardiovascular surgery, and reconstructive and aesthetic plastic surgery. Our instruments feature delicate design, industry-leading durability and a unique feel, which has made STILLE a well-known and highly respected supplier by surgeons throughout the world.

Surgical imaging tables

In our Surgical imaging tables business area, we offer specialty tables for minimally invasive procedures, with pricing, safety and features optimized for both in- and outpatient settings. The STILLE imagiQ2 surgical imaging table has the highest translucency in the world, and is optimized for advanced mobile Hybrid OR settings as well as peripheral and interventional procedures, enhancing crystal imaging while reducing radiation exposure as well as increasing efficiency in the OR.

The STILLE Medstone line of surgical imaging tables range from floating to fully motorized tables for vascular, electrophysiology, gastro (ERCP), pain management and urology to suit your needs, with reliability and optimal performance at an affordable cost at the highest standards of Swedish technology.



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30
Torshälla, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
sales@stille.se

www.stille.se

STILLE Surgical Inc.

915 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614
sales.us@stille.se