


STILLE


Surgical Instruments



Surgical perfection. For life.

Selected Models





STILLE Surgical Instruments Selected Models

Art. no: 950-G100-0001-EN-01.

Photos: Åke Gunnarsson, Vue and Jonas Frid, Studio Roland, Fanny Eriksson Palm

Printing: Arkitektkopia AB

© STILLE AB 2024

The world's best instruments, for the world's best surgeons

For more than 180 years, we have developed and manufactured the best surgical instruments for the world's most demanding surgeons. We would like to extend a sincere thank you to all our current customers and welcome our future customers to experience the STILLE difference. In this catalog we present a selection of STILLE original surgical instruments.

Precision, durability and feel are characteristic qualities of all STILLE instruments. The majority are handcrafted by our highly skilled artisans in Eskilstuna, Sweden. Our instrument manufacturing process from melting steel rods to finished product is a long one, involving multiple stages of development. Our method of crafting STILLE instruments involves only the highest-grade steel, which provides a unique feel and reliability to the surgeon in any surgical specialty.

In this catalogue you will find a wide selection of instruments, including the original SuperCut scissor developed in 1984 and our Diamond SuperCut design. Quality goes into our STILLE forceps, needle holders, clamps, rongeurs, curettes and retractors, as well as, a fine selection of stainless and titanium micro instruments. Common to all of our instruments is that they are smooth and consistently perform to your expectations. Our outstanding balance, solid design and ergonomics provides the surgeon with a reliable instrument throughout long and complex procedures.

We at STILLE are proud to have developed the world's first double-action bone rongeur. We also developed the world's first osteotome with the classic 8-sided handle. Many STILLE instruments have been copied, but a copy cannot measure up to the original.



Dr. Mark L. Jewell

My patients deserve nothing but the best. STILLE has been the natural extension of my hands for over 30 years.



– Dr. Mark L. Jewell, M.D.
Plastic Surgeon
www.markjewellmd.com

Table of Contents

6	Scissors	20	Clamps
10	Forceps	24	Skin Hooks and Retractors
13	Artery Forceps	25	Bone Instruments
14	Needle Holders	28	Micro Instruments
17	Bone Instruments		

SYMBOLS



SuperCut



Diamond Jaws



Ergonomic



TC-inserts



STILLE EYE SCISSORS

Art. Number: 101-8003

Delicate eye scissors with sharp or blunt tips and large finger rings. These are also called iris scissors and are used for nerves, tendons and vessels.



STILLE SUPERCUT EYE SCISSORS

Art. Number: 101-8015-1

SuperCut eye scissors with sharp or blunt tips and extremely sharp edges. Beveled blades.



KAYE STILLE SUPERCUT DISSECTING SCISSORS

Art. Number: 101-8021-1

Scissors with serrated blades for "holding" the tissue, making it easier to cut exactly as far as intended. SuperCut version.



JABALEY STILLE SUPERCUT ERGONOMIC DISSECTING SCISSORS

Art. Number: 101-8060-2

SuperCut scissors for dissection of tough tissue in hand surgery. Ergonomic model to avoid fatigue in the hand. Beveled blades.



MAYO STILLE DISSECTING SCISSORS

Art. Number: 101-8080-17

The most used scissors, for all types of surgery. Used for routine dissection and cutting in tissue. Beveled blades.



BABY METZENBAUM STILLE SUPERCUT DISSECTING SCISSORS

Art. Number: 101-8163

SuperCut dissecting scissors for tough tissue. Used in delicate surgery requiring small instruments.



METZENBAUM STILLE SUPERCUT DISSECTING SCISSORS

Art. Number: 101-8175-18

The popular Metzenbaum dissecting scissors with SuperCut edges. Rounded/beveled blades.



METZENBAUM STILLE DELICATE DIAMOND SUPERCUT DISSECTING SCISSORS

Art. Number: 101-8173-15-C

Delicate version with Super Cut edges and fine rounded tips.



**RAGNELL STILLE SUPERCUT
DISSECTING SCISSORS**

Art. Number: 101-8380-151

Dissecting scissors with flat, rounded tips for precise undermining of flaps. Often used in plastic surgery. SuperCut version.



**STEVENS-STILLE SUPERCUT
TENOTOMY SCISSORS**

Art. Number: 101-8167-151

Very popular SuperCut scissors for precise cutting and dissection of fine tissues. Rounded blades.



**STILLE SUPERCUT VASCULAR
SCISSORS**

Art. Number: 101-8200-60

SuperCut scissors with angled blades for better reach of vessels etc. Used in vascular surgery.



**PECK-JOSEPH STILLE NASAL
SCISSORS**

Art. Number: 101-8125

Double-beveled scissors for rhinoplasty.



**MATARASSO STILLE SUPERCUT
FACE-LIFT SCISSORS**

Art. Number: 101-8458-20

SuperCut face-lift scissors with double-beveled semi-sharp outer blades.



**LITTLER STILLE SUTURE
CARRYING SCISSORS**

Art. Number: 101-8164

Suture carrying scissors. Rounded blades.

ADSON STILLE TISSUE FORCEPS

Art. Number: 102-9022

Tissue forceps with very delicate points and 1x2 fine teeth. Used in most specialties for delicate work.

**GILLIES STILLE TISSUE FORCEPS**

Art. Number: 102-9027

Tissue forceps with 1x2 fine teeth and stabilizing pin. An excellent forceps for the surgeon who wants an Adson-size tip, but with greater length.

**STILLE 2X3 TEETH TISSUE FORCEPS**

Art. Number: 102-9062-20

Tissue forceps with 2x3 teeth. Used when more grasping power than 1x2 teeth is needed.



WAUGH-BROPHY TISSUE FORCEPS

Art. Number: 102-9252-25

Tissue forceps with 1x2 fine teeth. Very fine tips, like Adson. Useful when long and slender instruments are needed.



MCINDOE STILLE DRESSING FORCEPS

Art. Number: 102-9025-1

Dressing forceps with stabilizing pin and delicate jaws for grasping fine tissue.



STILLE SQUARE POINTS DRESSING FORCEPS

Art. Number: 102-9030-15

Dressing forceps with square shaped and serrated points used in all surgical specialties.



CRAFOORD-STILLE DRESSING FORCEPS

Art. Number: 102-9077-20

Dressing forceps with rounded points. The perfect balance and delicate grip gives this instrument a unique feeling.

**OLIVECRONA-STILLE DURA FORCEPS**

Art. Number: 102-9244-15

Dura forceps with 2x3 fine teeth. The nonpenetrating tips give the surgeon maximum grip without damaging the tissue.

**STILLE DE BAKEY VASCULAR FORCEPS**

Art. Number: 102-9312

Atraumatic forceps with De Bakey jaws for handling of blood vessels and other delicate tissue.

**BJÖRK-STILLE INSULATED DIATHERMY FORCEPS**

Art. Number: 102-9845-20

Insulated forceps for monopolar diathermy. Serrations allow grasping of delicate structures.





PEAN STILLE ARTERY FORCEPS

Art. Number: 103-10015-131

General utility forceps used in many surgical disciplines.



CRILE STILLE ARTERY FORCEPS

Art. Number: 103-10182-161


General purpose hemostat without teeth. Used in relatively small operating fields, when more strength is needed than a Mosquito can offer.



HALSTED-MOSQUITO STILLE ARTERY FORCEPS

Art. Number: 103-10032-121

Hemostat without teeth, used in small operative fields. For plastic surgery, neuro, hand or any delicate surgery.



LERICHE STILLE ARTERY FORCEPS

Art. Number: 103-10192

Fairly tapered hemostat. For use in small operative fields.



TC

DERF STILLE NEEDLE HOLDER*Art. Number: 104-15542*

Needle holder with rather short TC jaws, popular in plastic surgery. Intended for 3/0 suture and finer.



TC

WEBSTER STILLE NEEDLE HOLDER*Art. Number: 104-15559*

Delicate needle holder for fine needles. With smooth jaws or TC jaws. The Webster-Kleinert (104-15557) has offset bows for a comfortable grip. Intended for 4/0 suture and finer.



TC

HALSEY STILLE NEEDLE HOLDER*Art. Number: 104-15560*

Needle holder with TC jaws used in many procedures. One model with ergonomically designed finger ring for thumb. For slightly larger needles than the Webster needle holder. Intended for 3/0 suture and finer.



TC

MAYO-HEGAR STILLE NEEDLE HOLDER*Art. Number: 104-15561-181*

Commonly used general purpose needle holder with TC jaws. Intended for 0 suture and finer.



TC

CRILE-WOOD STILLE NEEDLE HOLDER*Art. Number: 104-15565-201*

General purpose needle holder with TC jaws. Intended for 2/0 suture and finer.



TC

CRILE-WOOD STILLE DELICATE NEEDLE HOLDER*Art. Number: 104-15566-151*

Delicate version of the popular Crile-Wood needle holder. With TC jaws. Intended for 3/0 suture and finer.



TC

OLSEN-HEGAR STILLE NEEDLE HOLDER*Art. Number: 104-15590-131*

Popular needle holder with scissors for suturing and cutting with the same instrument. Intended for 0 suture and finer. 3+3 ratchets, 104-15590-171 has 2+2 ratchets.



TC

STILLE-FRIMAND DUO-GRIP NEEDLE HOLDER*Art. Number: 104-15665-151*

Designed in cooperation with Carl-Fredrik Frimand, M.D. Intended for precise use in the palmar grip without the need of grip change during suturing. This new and innovative needle holder reduces the suturing time and tissue manipulation as well as offering better ergonomics for the surgeon. Intended for 2/0 suture and finer.

**RAPP-STILLE KNOT-HOLDING FORCEPS***Art. Number: 104-15140*

A 13 cm knot-holding forceps with smooth jaws. Ideal instrument to use as a tool for holding the suture knot in narrow areas.



STILLE EXTRACTION PLIERS

Art. Number: 105-16020

Original STILLE double-action joint facilitates removal of pins and wires.



STILLE BONE CUTTING FORCEPS CURVED

Art. Number: 105-16811

Original STILLE double-action angled bone cutting forceps.



BEYER STILLE BONE RONGEUR

Art. Number: 105-17210

Bone rongeur with original STILLE double-action joint and two different jaw widths.



LEKSELL-STILLE BONE RONGEUR

Art. Number: 105-17260

Bone rongeur with original STILLE double-action joint and sideways slightly curved jaws.



SEMB STILLE BONE HOLDING FORCEPS

Art. Number: 105-16150

Bone holding forceps with or without locking mechanism.



STILLE DOUBLE-ACTION JOINT WIRE CUTTER

Art. Number: 105-16223

Double joint wire cutter with TC-inserts. For hard wire up to 1.7 mm and soft wire up to 2.0 mm.



KERRISON STILLE LAMINECTOMY PUNCH

Art. Number: 105-17643-2-C

Kerrison punches with a unique rubber coated handle for improved ergonomics. Also available with coating for increased durability.



STILLE MULTI-PURPOSE CLAMP

Art. Number: 107-10825

Multi purpose clamps with angled, atraumatic DeBakey jaws.



SATINSKY STILLE TANGENTIAL OCCLUSION CLAMP

Art. Number: 107-10841

Tangential occlusion clamps with atraumatic DeBakey jaws.



STILLE DEBAKEY AORTIC-ANEURYSM CLAMP

Art. Number: 107-10849

Aortic aneurysm clamp with atraumatic DeBakey jaws.



STILLE ANASTOMOSIS CLAMP

Art. Number: 107-10860

Anastomosis clamps with atraumatic DeBakey jaws.



ERICSSON-STILLE CAROTID CLAMP

Art. Number: 107-10880

Carotid clamp with angled atraumatic DeBakey jaws.



BACKHAUS STILLE TOWEL CLAMPS

Art. Number: 107-10910

Towel clamps in two different versions. Sharp/sharp or ball/socket for paper towels.



BABCOCK STILLE INTESTINAL GRASPING FORCEPS

Art. Number: 107-12370-16

Intestinal grasping forceps in two different lengths.



**ALLIS STILLE INTESTINAL GRASPING
FORCEPS**

Art. Number: 107-13120-15

Tissue-grasping forceps for secure grip of tissue.



SCHROEDER STILLE TENACULUM FORCEPS

Art. Number: 107-13280

Uterine tenaculum forceps.

GILLIES STILLE SKIN HOOK

Art. Number: 108-24110

Skin hook with small or large hook on long flat shaft.



JOSEPH STILLE SKIN HOOK

Art. Number: 108-25012-10

Skin hook with sharp prongs. Available in five different models. Popular in plastic surgery.



STILLE BREAST RETRACTOR

Art. Number: 108-26038

Useful double-sided retractor for breast augmentation. Serrated blades in two different sizes and shapes. Ergonomic handle for maximum comfort.



WEITLANER STILLE RETRACTOR

Art. Number: 108-26254

General purpose self-retaining retractor used in many procedures.



STILLE BONE CURETTE

Art. Number: 111-20502-4

Oval and sharp bone curettes with widths from 2 mm to 12 mm.

**WIBERG STILLE RASPATORY**

Art. Number: 111-22171-7

A straight raspatory in four different sizes.

**HOWARTH STILLE NASAL RASPATORY**

Art. Number: 111-22370

Double-ended nasal raspatory with sharp hooks.

**STILLE NASAL RASPATORY**

Art. Number: 111-22520-4

Sharp, slightly curved nasal raspatory.



FREER STILLE ELEVATOR

Art. Number: 111-22140-6

Double-ended elevator sharp/blunt with rounded, delicate handle.



ORIGINAL STILLE DELICATE OSTEOTOME

Art. Number: 111-21305-8

The famous STILLE original osteotome, available in 5 different sizes.



STILLE HAND OSTEOTOME

Art. Number: 111-21456-4

Delicate osteotome with phenolic handle. Intended for hand surgery, but also useful in other fine orthopaedic work.



LEXER STILLE OSTEOTOME

Art. Number: 111-21460-15

Osteotome with phenolic handle.

**ORIGINAL STILLE GOUGE**

Art. Number: 111-21315-8

Original STILLE gouge. Available in five different sizes from 6 mm to 15 mm.

**STILLE STAINLESS STEEL MALLET**

Art. Number: 111-21820

Stainless steel mallet.





STILLE SUPERCUT MICRO SCISSORS

Art. Number: 117-A110-18

Fine micro scissors with straight or curved blades with sharp or blunt tips. All models with SuperCut blades.



STILLE SUPERCUT MICRO VASCULAR SCISSORS

Art. Number: 117-A120-45

SuperCut micro scissors available in five different blade angles for fine vascular surgery. Blade length 11 mm.



STILLE MICRO FORCEPS

Art. Number: 117-B110-18

Well balanced and ergonomic micro forceps of ultra high quality. With tying platform and straight or curved jaws. 0.4 mm tips.





STILLE MICRO RING FORCEPS

Art. Number: 117-B140-18

Counter-balanced micro forceps with fine Diamond Jaws ring tips.



STILLE MICRO DEBAKEY FORCEPS

Art. Number: 117-B160-18

Round handle, with 1.2 mm fine atraumatic DeBakey tips.



STILLE MICRO NEEDLE HOLDER REGULAR BOX LOCK WITH RATCHET

Art. Number: 117-C121-18

Round handle, fine 0.8 mm Diamond Jaws and Regular Box Lock for sutures 7-0 and smaller. With ratchet.



STILLE TITANIUM MICRO FORCEPS

Art. Number: 117-B120-18-TI

Round handle, with tying platform, straight or curved 0.4 mm tips.



**STILLE TITANIUM COUNTER-BALANCED MICRO FORCEPS**

Art. Number: 117-B130-21-TI

Round handle, with tying platform, counter-balanced, straight or curved 0.4 mm tips.

**STILLE TITANIUM MICRO RING FORCEPS**


Art. Number: 117-B140-21-TI

Round handle, with fine ring tips, counter-balanced.

**STILLE TITANIUM MICRO NEEDLE HOLDER REGULAR BOX LOCK**

Art. Number: 117-C110-18-TI


Round handle, fine 0.8 mm Diamond Jaws and Regular Box Lock for sutures 7-0 and smaller. With or without ratchet.

**STILLE TITANIUM MICRO NEEDLE HOLDER STREAMLINE BOX LOCK**

Art. Number: 117-C125-21-TI

Round handle, 1.2 mm Diamond Jaws and Streamline Box Lock for sutures 5-0 and smaller. With or without ratchet.



The background is a dark blue gradient. Overlaid on it are white line-art illustrations of various surgical instruments. In the upper right, there is a pair of large surgical scissors with long handles and curved blades. On the left side, there is a pair of long-handled surgical forceps with multiple circular joints. At the bottom, there is a pair of long, thin surgical forceps with a textured grip section on the handle.

**“An investment in a STILLE
instrument is an investment that
will last throughout your career.”**

– Dr. Jordan Deschamps-Braly, M.D.
Plastic Surgeon

Surgical perfection. For life.

**For a complete presentation
of all STILLE instruments,**

Scan the QR code to visit our website and
see our product catalogue



950-G100-0001-EN-01

STILLE AB

Ekbacken 11
SE-644 30 Torshälla, Sweden
☎ +46 8 588 580 00
✉ info@stille.se
🌐 www.stille.se

STILLE SURGICAL INC.

840 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614



Surgical perfection. For life.