




Surgical perfection. For life.



# Are your instruments safe to use?



Surgical instruments are affected by age, usage and handling. This means that many of the instruments on your trays can be dangerous to use, as there may be defects that are difficult to see, e.g. corrosion and small fractures. We can help you to find those instruments.

STILLE has manufactured surgical instruments of the absolutely best quality since 1841. Still today, just like 175 years ago, the instruments are mainly manufactured by hand, all the way from forging to shaping, grinding, polishing and adjustment. Probably nobody knows better than STILLE what characterizes an instrument of good quality and endurance, or how correct maintenance can increase its life expectancy and contribute to better patient safety.

Correct maintenance and care drastically increases the instrument's life expectancy. By taking good care of the instruments, the clinic can use their precious instruments

longer, transforming the investments in expensive instruments to better long term economy.



We categorize the instruments in three groups;

- 1. Instruments that are in good shape and safe to use without any remedy.**
- 2. Instruments that are in need for service. These instruments are marked to enable sending for service following a schedule that matches the clinic's operations.**
- 3. Instruments that are not safe to use and need to be replaced. These are specifically marked and put aside.**

Following to the instrument audit, you will receive a written report with exact numbers of each category as well as recommendations how to treat the instruments to avoid problems as we can see during the audit.

We can assist you with a professional instrument audit, made by one of our skilled instrument makers with many years of experience. You are welcome to contact our local STILLE-partner to learn more and find a suitable setup and date for your audit. Calculate with appr 1.200 instruments per day.

## STILLE AB

Ekbacken 11  
SE 644 30 Torshälla  
Telephone: +46 8 588 580 00  
sales@stille.se

[www.stille.se](http://www.stille.se)