



Stille *Visi*™ XLT-125

Freedom of movement



Surgical Perfection. For life.

Stille Visi™

It just brings perfect light

Stille Visi headlight provides the perfect light source and offers full mobility for gynecology, urology and general surgery procedures such as orthopedic, plastic, neuro and cardiovascular surgery. Stille Visi XLT-125 has adjustable brightness ranging from 40 000 all the way up to 125 000 Lux from one single LED. Adjustable aperture ensures the perfect light for every procedure. Stille provides a 5-year LED warranty.

Minimal maintenance. There are no cables to change and the LED will last for at least 50 000 hours without losing its brightness.

Flexible arm allows easy use with or without loupes.

White light (6 100 K) for optimal display of a range of tissues.

Adjustable aperture 5-15 cm on normal working distance.

Ergonomic and light weight headlight part of only 198 g.

StilleVisi XLT-125
Headlight System
215-90071

Rechargeable B2X battery pack with a battery life of 4-12 hours of continuous surgery depending on chosen light intensity. A B4X battery pack is offered as accessory which offers a battery life of 8-24 hours!

Accessories

- 215-80077-001 B4X, Battery Unit, big, with clips function
- 215-80076-001 B2X, Battery Unit, small, with clips function
- 215-80033-002 Set of Comfort Pads
- 215-61507-001 Belt for Battery Unit
- 215-80031-002 Hard Plastic Carrying Case
- 215-80031-003 Soft Plastic Carrying Case
- 215-80058-001 Battery Charger with Cable
- 215-80042-001 Battery Cable



215-80076-001
B2X, Battery Unit



215-80077-001
B4X, Battery Unit

StilleVisi XLT-125 Headlight System

Includes: 1 Headlight adjustable 40 000 -125 000 lux, B2X Battery Unit with clips function, Battery Unit Charger, Belt, 1 Soft Plastic Carrying Case.

Article code 215-90071



Surgical Perfection. For life.

Stille AB

Gustav III:s Boulevard 42
SE-169 73 Solna, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
sales@stille.se

Distributed by: