

High Quality Hospital Grade Surgical Imaging Table with Free Float Capabilities

The STILLE Medstone surgical imaging tables contribute to enhancing efficiency in the OR. Representing value priced solutions with quality and reliability.



Standard configuration and specifications

- Tabletop dimensions 24–82"
- Heavy duty steel base frame
- Four hospital grade locking casters
- Float control unit
- Mattress pad 2" thick with "Softcare™" technology
- One (1) hand control
- One (1) foot control
- One (1) catheter tray 24"
- Three (3) Patient safety straps
- 110VAC, 230VAC available

Optional Accessories

- Adjustable patient arm board with pad (shown)
- Fistula arm board with pad
- Radiation protection shield
- Short side rail for accessories

The STILLE Medstone3 CV table with free floating capabilities allows for head to toe imaging with complete accuracy, speed and freedom of movements. Ideal for endovascular and fluoroscopy imaging.

KEY features

- Longitudinal float 36"
- Transverse float ± 4 "
- Imaging length 64" head to toe
- 3D imaging length 32"
- Motorized elevation 26–42"
- Patient weight capacity 500lbs
- Parts and labor warranty 3/2 year



Optional basic arm board

STILLE Medstone Free Float Imaging Table

Part number	Model reference	Description
500-1025	Medstone3 CV	Vascular table 36" longitudinal, transverse $\pm 4"$, 26-42" elevation, imaging length 64" head to toe
505-1116	Optional	Radiation protection shield
505-1117	Optional	Basic arm board (mounted directly onto tabletop)
505-1118	Optional	Fistula arm board
505-1119	Optional	Short rail 6"
505-1121	Optional	Pair of clamps to attach catheter tray to head end of table

Specification tolerance: $\pm 1^\circ$, $\pm .5"$.



STILLE Medstone | Vascular

STILLE – Surgical perfection since 1841

STILLE is one of the oldest medical device companies in the world founded at Karolinska University Hospital in Sweden back in 1841. At STILLE, we develop, produce and sell high-quality surgical instruments and surgical imaging tables to physicians worldwide.

STILLE's offering consists of a strong brand and products renowned for their superior quality. We manage our business in two distinct business areas—Surgical Instruments and Surgical imaging tables.

Surgical instruments

In our Surgical Instruments business area, we have a complete offering of handcrafted, forged instrumentation for all fields of open surgery, with a special emphasis on cardiovascular surgery, and reconstructive and aesthetic plastic surgery. Our instruments feature delicate design, industry-leading durability and a unique feel, which has made STILLE a well-known and highly respected supplier by surgeons throughout the world.

Surgical imaging tables

In our Surgical imaging tables business area, we offer specialty tables for minimally invasive procedures, with pricing, safety and features optimized for both in- and outpatient settings. The STILLE imagiQ2 surgical imaging table has the highest translucency in the world, and is optimized for advanced mobile Hybrid OR settings as well as peripheral and interventional procedures, enhancing crystal imaging while reducing radiation exposure as well as increasing efficiency in the OR.

The STILLE Medstone line of surgical imaging tables range from floating to fully motorized tables for vascular, electrophysiology, gastro (ERCP), pain management and urology to suit your needs, with reliability and optimal performance at an affordable cost at the highest standards of Swedish technology.



Surgical Perfection

STILLE AB

Ekbacken 11
SE-644 30
Torshälla, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
sales@stille.se

STILLE Surgical Inc.

915 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614
sales.us@stille.se

www.stille.se