



Surgical perfection. For life.

STILLE Medstone ERCP / Endoscopic fluoroscopy

Standard Warranty

3 years parts / 3 years labor

High Quality Hospital Grade Surgical Imaging Tables with Extra Wide Table Top for ERCP and other Endoscopic Procedures

The STILLE Medstone C-arm tables contribute to enhancing efficiency while representing an affordable solution with quality and reliability provided by Swedish technology.

STILLE Medstone ERCP Surgical Imaging Tables

3- and 5- movement C-arm tables with extra wide table top.

- Rectangular 26x82" carbon fiber table top without facial cutout
- Weight capacity 500 lbs
- Motorized movements
- Four hospital grade locking castors
- Metal free imaging 58"
- Carbon fiber table top with <1 mm Aluminium Equivalence
- Parts and labor warranty 3 years
- 110VAC, 60Hz
- 2" table mattress and 3 patient straps
- Wired hand control

3 and 5 Motorized Movement Options Available:

- Elevation
- Lateral tilt
- Trendelenburg
- Longitudinal travel (Only Medstone5 ERCP)
- Transverse travel (Only Medstone5 ERCP)



Photo shows Medstone5 ERCP (500-1031) with infectious control covers

Model	Medstone3 ERCP	Medstone5 ERCP
Part number	500-1040*	500-1031*
Table top	26" x 82"	26" x 82"
Elevation Range	26.6"-40.8"	26.3"-40.1"
Lateral tilt	+/- 15°	+/- 15°
Trendelenburg	+/- 15°	+/- 15°
Longitudinal travel	x	15.75"
Transverse travel	x	6"
Infectious Control Covers	x	Yes

Specification tolerance: $\pm 1^\circ$, $\pm .5"$

*Basic configuration without detachable side rails

Optional Accessories

Part number	Description
505-1107	Three (3) movement foot control for elevation, Trendelenburg and lateral tilt (fits table 500-1040)
505-1109	Five (5) movement foot control for tables with elevation, Trendelenburg, lateral tilt, longitudinal travel, transverse travel (fits table 500-1031)
535-1770	X-ray protective kit grey, 6 single panels; attachment to table top with self-adhesive Velcro © strips
505-1137	Extension cable hand control

STILLE – Surgical perfection. For life.

STILLE is one of the oldest medical device companies in the world. At STILLE, we develop, produce and sell high-quality surgical instruments and surgical imaging tables to physicians worldwide, and has done so for 170 years.

STILLE's offering consists of a strong brand and products renowned for their superior quality. We manage our business in two distinct business areas – Surgical Instruments and Surgical imaging tables.

Surgical instruments

In our Surgical Instruments business area, we have a complete offering of handcrafted, forged instrumentation for all fields of open surgery, with a special emphasis on cardiovascular surgery, and reconstructive and aesthetic plastic surgery.

Our instruments feature delicate design, industry-leading durability and a unique feel, which has made STILLE a well-known and highly respected supplier by surgeons throughout the world.

Surgical imaging tables

In our Surgical imaging tables business area, we offer specialty tables for minimally invasive procedures, with pricing, safety and features optimized for both in- and outpatient settings. The STILLE imagiQ2 surgical imaging table has the highest translucency in the world, and is optimized for advanced mobile Hybrid OR settings as well as peripheral and interventional procedures.

The STILLE Medstone line of imaging tables offer reliability and optimal performance at an affordable cost with the highest standards of Swedish technology.



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30
Torshälla, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
sales@stille.se

www.stille.se/us

STILLE Surgical Inc.

915 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614
sales.us@stille.se