



Stille Ergo Scissors

More comfort, same precision



Surgical Perfection. For life.

Reducing fatigue, enhancing comfort



Ergonomic design that maintains your performance during time-consuming procedures, without compromising precision.

Surgical scissors are normally held with precision grip using the thumb and the long or ring finger. When using scissors of a traditional design, the hand usually ends up in a strenuous position which may cause tension or fatigue in the surgeon's hand after sustained and repeated use.

In order to achieve a grip that more closely resembles the hand's normal resting position, Stille has developed a range of scissors with increased space between the finger rings compared to traditional scissors. Stille's Ergo scissors have the advantage to avoid fatigue in the surgeon's hand during time-consuming procedures. As a result, the hand joint is painless and the risk for arthritis is minimized.

The dissection work becomes easier and with less effort.

For more information please visit: www.stille.se



Stille Ergo product range

Description	Article code
STILLE BLEPHAROPLASTY ERGO SUPER CUT SCISSORS curved 9,5 cm 3,7 5 in	101-8019-2
STILLE RAGNELL ERGO SUPER CUT SCISSORS curved 9.5 cm 3.75 in	101-8031-2
STILLE JABALEY ERGO SUPER CUT SCISSORS straight 12.5 cm 5 in	101-8060-2
STILLE JABALEY ERGO SUPER CUT SCISSORS curved 12.5 cm 5 in	101-8062-2
STILLE MAYO ERGO SUPERCUT SCISSORS rounded blades curved 15 cm 6 in	101-8087-152
STILLE MAYO ERGO SUPERCUT SCISSORS rounded blades curved 17 cm 6.75 in	101-8087-172
STILLE BABY METZENBAUM ERGO SUPER CUT SCISSORS curved 12.5 cm 5 in	101-8163-2
STILLE STEVENS TENOTOMY ERGO SUPER CUT SCISSORS straight 12.5 cm 5 in	101-8166-2
STILLE STEVENS TENOTOMY ERGO SUPER CUT SCISSORS curved 12.5 cm 5 in	101-8167-2
STILLE STEVENS TENOTOMY ERGO SUPER CUT SCISSORS curved 15 cm 6 in	101-8167-152
STILLE STEVENS TENOTOMY ERGO SUPER CUT SCISSORS curved 18 cm 7 in	101-8167-182
STILLE METZENBAUM ERGO SUPER CUT curved 15 cm 6 in	101-8175-152
STILLE METZENBAUM ERGO SUPER CUT curved 18 cm 7 in	101-8175-182
STILLE RAGNELL ERGO SUPERCUT curved 12.5 cm 5 in	101-8380-122
STILLE RAGNELL ERGO SUPERCUT curved 15 cm 6 in	101-8380-152
STILLE GORNEY FACELIFT ERGO SCISSORS SUPER CUT straight 17 cm 6.75 in	101-8430-172
STILLE GORNEY FACELIFT ERGO SCISSORS SUPER CUT curved 17 cm 6.75 in	101-8431-172
STILLE GORNEY FACELIFT ERGO SCISSORS SUPER CUT curved 20 cm 8 in	101-8433-202
STILLE GORNEY FACELIFT ERGO SCISSORS SUPER CUT serrated curved 20 cm 8 in	101-8435-202
STILLE-FREEMAN FACELIFT ERGO SCISSORS SUPER CUT curved 18 cm 7 in	101-8445-18
STILLE MATARASSO FACELIFT ERGO SCISSORS SUPERCUT curved 17 cm 6.75 in	101-8458-172
STILLE MATARASSO FACELIFT ERGO SCISSORS SUPERCUT curved 20 cm 8 in	101-8458-202



Stille Ergo scissors

- ✓ Ergonomically designed to avoid fatigue
- ✓ 30 year "career" warranty on all Stille produced instruments
- ✓ Free service voucher included





| **Stille Ergo Scissors** |

Stille - Surgical Perfection

Stille is one of the oldest medical device companies in the world. We develop, manufacture and sell proprietary high-quality surgical instruments and specialty tables to leading surgeons worldwide. We have done so for more than 170 years.

Stille's offering consists of a strong brand and products renowned for their superior quality. We manage our business in two distinct business areas—Surgery and Patient Positioning.

Stille surgical instruments are considered state of the art by many surgeons around the world. All Stille instruments are handcrafted as an art by highly skilled instrument makers in Eskilstuna, Sweden.

The slender design and excellent balance are the characteristics which make the Stille Instruments easily recognizable. Their highly polished surface and delicate design make Stille instruments lighter and easier to keep clean, thus more durable against discoloration and corrosion. All this provides Stille instruments with a unique and outstanding feeling.



Surgical Perfection. For life.

Stille AB

Gustav III:s Boulevard 42
SE-169 73 Solna, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
info@stille.se

www.stille.se

Stille Surgical Inc.

943 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614

Distributed by: