



Nålförare storlekstabell | Needle holder in relation to suture size | Nadelhalter in Relation zur Stärke des Nahtmaterials |
 Il portaaghi in relazione alla misura della sutura | Porta-agujas en relación al tamaño de la sutura | 縫合糸サイズに適する持針器



Art.no	Description	Max suture size
104-15542	RAPP-STILLE KNOT-HOLDING FORCEPS 13 cm 5"	3/0
104-15545	WEBSTER STILLE very delicate DD jaws 12 cm 4.75"	5/0
104-15546	WEBSTER STILLE very delicate smooth 12 cm 4.75"	5/0
104-15549-131	KILNER STILLE TC jaws 13 cm 5"	4/0
104-15550	GILLIES STILLE serrated 16 cm 6.25"	2/0
104-15551	GILLIES STILLE TC jaws 16 cm 6.25"	2/0
104-15556	WEBSTER STILLE delicate smooth 12 cm 4.75"	4/0
104-15557	WEBSTER-KLEINERT STILLE smooth offset bows 12 cm 4.75"	4/0
104-15558	WEBSTER STILLE smooth TC jaws 12 cm 4.75"	3/0
104-15559	WEBSTER STILLE serrated TC jaws 12 cm 4.75"	4/0
104-15560	HALSEY STILLE TC jaws 12.5 cm 5"	3/0
104-15560-2	HALSEY STILLE ERGONOMIC TC jaws 12.5 cm 5"	3/0
104-15561-141	MAYO-HEGAR STILLE TC jaws 14 cm 5.5"	0
104-15561-161	MAYO-HEGAR STILLE TC jaws 16 cm 6.25"	0
104-15561-181	MAYO-HEGAR STILLE TC jaws 18 cm 7"	0
104-15561-201	MAYO-HEGAR STILLE TC jaws 20 cm 8"	0
104-15561-251	MAYO-HEGAR STILLE TC jaws 25 cm 10"	0
104-15562-161	MAYO-HEGAR STILLE Heavy Pattern TC jaws 16 cm 6.25"	1
104-15562-181	MAYO-HEGAR STILLE Heavy Pattern TC jaws 18 cm 7"	1
104-15562-201	MAYO-HEGAR STILLE Heavy Pattern TC jaws 20 cm 8"	1
104-15562-241	MAYO-HEGAR STILLE Heavy Pattern TC jaws 24 cm 9.5"	1
104-15562-271	MAYO-HEGAR STILLE Heavy Pattern TC jaws 27 cm 10.5"	1
104-15562-301	MAYO-HEGAR STILLE Heavy Pattern TC jaws 30 cm 12"	1
104-15565-151	CRILE-WOOD STILLE TC jaws 15 cm 6"	2/0
104-15565-181	CRILE-WOOD STILLE TC jaws 18 cm 7"	2/0
104-15565-201	CRILE-WOOD STILLE TC jaws 20 cm 8"	2/0
104-15565-231	CRILE-WOOD STILLE TC jaws 23 cm 9"	2/0

Art.no	Description	Max suture size
104-15566-151	BABY CRILE-WOOD STILLE delicate TC 15 cm 6"	3/0
104-15566-181	CRILE-WOOD STILLE delicate TC jaws 18 cm 7"	3/0
104-15566-201	CRILE-WOOD STILLE delicate TC jaws 20 cm 8"	3/0
104-15566-231	CRILE-WOOD STILLE delicate TC jaws 23 cm 9"	3/0
104-15571	NEW ORLEANS STILLE TC jaws 18 cm 7"	3/0
104-15573-161	RYDER STILLE TC jaws 16 cm 6.25"	3/0
104-15573-181	RYDER STILLE TC jaws 18 cm 7"	3/0
104-15573-201	RYDER STILLE TC jaws 20 cm 8"	3/0
104-15585-161	STILLE DD jaws 16 cm 6.25"	4/0
104-15585-181	STILLE DD jaws 18.5 cm 7.5"	4/0
104-15585-221	STILLE DD jaws 22 cm 8.75"	4/0
104-15585-251	STILLE DD jaws 25 cm 10"	4/0
104-15586	STILLE TC jaws 18 cm 7"	4/0
104-15590-121	OLSEN-HEGAR STILLE TC jaws with scissors 12 cm 4.75"	0
104-15590-131	OLSEN-HEGAR STILLE TC jaws with scissors 13.5 cm 5.25"	0
104-15590-171	OLSEN-HEGAR STILLE TC jaws with scissors 17 cm 6.75"	0
104-15590-181	OLSEN-HEGAR STILLE TC jaws with scissors 18.5 cm 7.25"	0
104-15621-179	MATHIEU STILLE TC jaws 17 cm 6.75"	1
104-15621-209	MATHIEU STILLE TC jaws 18.5 cm 7.25"	1
104-15777	KAYE STILLE with scissors serrated TC jaws 12 cm 4.75"	5/0
104-15785	CONVERSE STILLE TC jaws 10.5 cm 4.25"	4/0

Suture size indicates the largest recommended for each needle holder.