



Surgical perfection. For life.

STILLE Medstone URO/ GYN Urology/ Cysto/ Brachytherapy/ Gynecology

High Quality Hospital Grade C-Arm Tables for Urology and Gynecology

The STILLE Medstone c-arm tables contribute to enhancing efficiency while representing an affordable solution with quality and reliability provided by Swedish technology.

STILLE Medstone URO Surgical Imaging Table

4- movement C-arm tables with carbon fiber table top and optional accessories.

- Rectangular 61 cm wide table top
- Weight capacity 227 kg
- Motorized movements
- Four hospital grade locking castors
- Imaging area 61x61 cm
- Carbon fiber table top with <1 mm Aluminium Equivalence
- 230Vac, 3,4A, 50 Hz and 120Vac, 6.5A, 60 Hz with 230Vac or 120Vac style line cord
- 5 cm table mattress and 3 patient straps
- Infectious control covers
- Wired hand control
- Wired foot control with pedestal to place the foot control on the table base
- One flexible uro-catcher drain bag mounting hoop
- One starter package of 20 disposable drain bags

4 Motorized Movements Available:

- Height adjustment
- Trendelenburg
- Longitudinal travel
- Lateral travel



Photo shows Medstone4 URO with leg holders 505-1126

Model	Medstone4 URO
Part number	500-1028
Table top	61x122 cm
Height adjustment	66-101 cm
Trendelenburg	± 15°
Lateral tilt	x
Longitudinal travel	40 cm
Lateral travel	15 cm
Auto horizontal and center home positioning	✓
Infectious Control Covers	✓

Specification tolerance: ±1°, ± 1.27 cm

Optional Accessories

Part number	Description
505-1126	Leg holders Goepel type/ pair with attachment clamps
505-1114	Leg holders (pair) with gas spring support and quick lock release handle
505-1125	Radiolucent leg extension plate 61x71 cm, 7.6 kg incl. pad
505-1170	Non-radiolucent leg extension plate 61x64 cm, 4.3 kg incl. pad
514-50-3SR025	Universal adjustable arm board with pad and clamp
505-1127	10 non-sterile drain bags
505-1120	Replacement uro-catcher drain bag mounting loop
535-1770	X-ray protective kit grey, 6 single panels; attachment to table top with self-adhesive Velcro © strips
535-1771	X-ray protective kit blue, 6 single panels; attachment to table top with self-adhesive Velcro © strips
535-1772	X-ray protective loop grey, 5 connected panels, attachment to side rails or table top with self-adhesive Velcro © strips
535-1773	X-ray protective loop blue, 5 connected panels, attachment to side rails or table top with self-adhesive Velcro © strips




505-1114
Leg holders with gas spring support (pair)



505-1126
Leg holders Goepel (pair)



514-50-3SR025
Adjustable arm board with pad



505-1125
Radiolucent leg extension plate



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30
Torshälla, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
sales@stille.se

www.stille.se