



STILLE Medstone ERCP

High Quality Hospital Grade Surgical Imaging Tables for Gastroenterology / ERCP

The STILLE Medstone C-arm tables contribute to enhancing efficiency while representing an affordable solution with quality and reliability provided by Swedish technology.

STILLE Medstone2 ERCP, 3 ERCP, 5 ERCP

- ▶ 2, 3 or 5 motorized movement options available
- ▶ Metal free imaging 147 cm / 58"
- ▶ Carbon fiber table top with <1 mm Al. Equivalence
- ▶ 66 x 208 cm / 26" x 81.9" rectangular table top
- ▶ 230 Vac, 3.4 A, 50 Hz with 230 Vac style line cord; certain models are available as a 120Vac US version, 6.5 A, 60 Hz with 120 Vac style line cord
- ▶ 5 cm / 2" table mattress and 3 patient straps
- ▶ Four hospital grade locking castors
- ▶ Wired hand control

- ✓ 227 kg / 500 lbs weight capacity
- ✓ With extra wide carbon fiber table top

	Medstone2 ERCP (500-1041)	Medstone3 ERCP (500-1007)	Medstone3 ERCP (500-1040)	Medstone3 ERCP (500-1044)	Medstone5 ERCP (500-1031)
Height adjustment⁽¹⁾ (cm)	✓ 66.5 - 101.9 (26.2" - 40.1")	✓ 63.5 - 98.8 (25.0" - 38.9")	✓ 67.7 - 103.7 (26.7" - 40.8")	✓ 63.9 - 98.9 (25.2" - 38.9")	✓ 67.0 - 102.0 (26.4" - 40.2")
Trendelenburg⁽²⁾ (± 15°)	✗	✗	✓	✓	✓
Lateral tilt⁽²⁾ (± 15°)	✓	✗	✓	✗	✓
Longitudinal travel⁽¹⁾ (40 cm)	✗	✓	✗	✓	✓
Lateral travel⁽¹⁾ (15 cm)	✗	✓	✗	✗	✓

⁽¹⁾ ± 1.25cm / 0.5" tolerance ; ⁽²⁾ ± 2.0° tolerance

Optional Accessories*



Medstone ERCP package (505-1189)

includes:

- ▶ 1 pair of 80.5 cm side rails (505-1181)
- ▶ 1 pair of side gates (514-100530)
- ▶ 1 adjustable arm board with mattress and clamp (514-50-3SR025)

Foot Control

- ▶ 505-1140: for table 500-1041
- ▶ 505-1139: for table 500-1007
- ▶ 505-1107: for table 500-1040
- ▶ 505-1138: for table 500-1044
- ▶ 505-1109: for table 500-1031



X-ray protective shields

- ▶ 535-1773: "loop blue". 5 connected panels; attachment to side rails or, to table top with self-adhesive Velcro® strips. Here attached to the side rail.
- ▶ 535-1771: "kit blue". 6 single panels; attachment to table top with self-adhesive Velcro® strips. Here attached to the head end of the table.

*Please see the accessory brochure for the full list (anesthesia screen, IV pole, ...).

STILLE AB
Ekbacken 11
644 30 Torshälla
Sweden

☎ +46 8 588 580 00
✉ sales@stille.se
🌐 www.stille.se

STILLE Surgical Inc.
840 Parkview Blvd
Lombard, IL 60148
USA

☎ +1 800 665 1614
✉ sales.us@stille.se
🌐 www.stille.se/us

These products may not be available in all regions. Please contact your STILLE-representative for more information on market availability.

This document correctly describe these devices at the time the document was printed. Please always refer to the user manual for any product amendments.

The product images shown are for illustration purposes only and may not be an exact representation of the product.