




Surgical perfection. For life.



NEW

STILLE imagiQ2™

3D Head Rest / Extension

Taking imaging to new dimensions.

STILLE imagiQ2™ 3D Head rest.

The latest accessory to enhance possibilities of imaging with a full carbon fiber extension.




Freedom of imaging.

- ✓ A must have for 3D rotations and obliques with 3D/2D C-arms
- ✓ Extends radiolucent length by 30cm / 11,8"
- ✓ Carbon Fiber laminated board with an antistatic mattress

Simple.

- ✓ Effortlessly slides on the head end of the table
- ✓ Easy to clean and disinfect
- ✓ Compatible only with short side rails
535-1724 (US) or 535-1725 (EU)



STILLE AB

Ekbacken 11
644 30 Torshälla
Sweden

☎ +46 8 588 580 00
✉ sales@stille.se
🌐 www.stille.se

STILLE Surgical Inc.

840 Parkview Blvd
Lombard, IL 60148
United States

☎ +1 800 665 1614
✉ sales.us@stille.se
🌐 www.stille.se/us/

This product may not be available in all regions. Please contact your STILLE-representative for more information on market availability. This document describes the device at the time it was printed. Please always refer to the user manual / instructions for use for more details and any product amendments. The product images shown are for illustration purposes only and may not be an exact representation of the product.