



Surgical perfection. For life.



STILLE

Popular Models: Needle Holders

A perfect suture. Every time.

SYMBOLS

Diamond Jaws



TC-inserts

STILLE offers the highest quality handmade instruments for the world's best surgeons

Few surgical instruments have to be more enduring than needle holders. STILLE's needle holders guarantee a firm and secure grip of the needle. They are easily recognized by their delicate and soft design. Once you try them, you will know the difference. STILLE enhances your investment in quality by providing you with a 30-year instrument warranty.

A perfect result when suturing requires a reliable and durable instrument. Stille needle holders are easily recognized by their delicate and soft design, which helps the plastic surgeon to hold the needle very securely. The golden finger loops indicate that the jaws are fitted with inserts of hard metal, which prolongs the life of the instrument. Some needle holders are instead coated with Diamond Jaws where hard metal is sprayed on as a thin layer on the inside of the jaws.

Another feature of Stille's needle holders is the open box lock. This reduces the risk of the needle holder being caught in the sutures or the surgeon's gloves. It also simplifies the cleaning and sterilization process. The ratchet's unique design ensures that the instrument stays closed in a safe position while remaining easy to open. All of these features are characteristic of that special feel that instruments from Stille provide.

For precise matching of sutures and needle holders, see the separate guide for matching suture size to individual needle holders.

KNOT HOLDING FORCEPS

The STILLE Rapp knot holding forceps is a well known and unique aid for holding the knot when tying the suture.



Needle Holders

POPULAR MODELS

- 6 Needle Holders
- 7 Knot Holding Forceps
- 9 Wire Twisting & Holding Forceps



104-15542

DERF STILLE NEEDLE HOLDER

Needle holder with rather short TC jaws, popular in plastic surgery. Intended for 3/0 suture and finer.

| ART. NUMBER | LENGTH | TC |
|-------------|-------------|----|
| 104-15542 | 12 cm 4.75" | • |



Art. Number: 104-15542



104-15545



104-15546

WEBSTER STILLE NEEDLE HOLDER

Very delicate needle holder for very fine sutures. Smooth jaws or Diamond jaws. Intended for 5/0 suture and finer.

| ART. NUMBER | JAWS | LENGTH | TC |
|-------------|----------------------|-------------|----|
| 104-15545 | Very delicate | 12 cm 4.75" | • |
| 104-15546 | Very delicate Smooth | 12 cm 4.75" | • |



Art. Number: 104-15545



104-15559

WEBSTER STILLE NEEDLE HOLDER

Delicate needle holder for fine needles with TC jaws. Intended for 4/0 suture and finer.

| ART. NUMBER | JAWS | LENGTH | TC |
|-------------|----------|-------------|----|
| 104-15559 | Serrated | 12 cm 4.75" | • |



Art. Number: 104-15559



Art. Number: 104-15560



104-15560

HALSEY STILLE NEEDLE HOLDER

Needle holder with TC jaws used in many procedures. For slightly larger needles than the Webster needle holder. Intended for 3/0 suture and finer.

| ART. NUMBER | LENGTH | TC |
|-------------|------------|----|
| 104-15560 | 12.5 cm 5" | • |



104-15561

MAYO-HEGAR STILLE NEEDLE HOLDER

Commonly used general purpose needle holder with TC jaws. Intended for 0 suture and finer.

| ART. NUMBER | LENGTH | TC |
|---------------|-------------|----|
| 104-15561-141 | 14 cm 5.5" | • |
| 104-15561-161 | 16 cm 6.25" | • |
| 104-15561-181 | 18 cm 7" | • |
| 104-15561-201 | 20 cm 8" | • |
| 104-15561-251 | 25 cm 10" | • |



Art. Number: 104-15561-181



Art. Number: 104-15562-201

MAYO-HEGAR STILLE HEAVY NEEDLE HOLDER

Commonly used general purpose needle holder with TC jaws. Intended for 0 suture and finer.

| ART. NUMBER | PATTERN | LENGTH | TC |
|---------------|---------|---------------|----|
| 104-15562-161 | Heavy | 16 cm 6.25" | • |
| 104-15562-181 | Heavy | 18.5 cm 7.25" | • |
| 104-15562-201 | Heavy | 20.5 cm 8" | • |



Art. Number: 104-15565-201



104-15565

CRILE-WOOD STILLE NEEDLE HOLDER

General purpose needle holder with TC jaws. Intended for 2/0 suture and finer.

| ART. NUMBER | LENGTH | TC |
|---------------|----------|----|
| 104-15565-151 | 15 cm 6" | • |
| 104-15565-181 | 18 cm 7" | • |
| 104-15565-201 | 20 cm 8" | • |
| 104-15565-231 | 23 cm 9" | • |



Art. Number: 104-15566-151



104-15566

CRILE-WOOD STILLE DELICATE NEEDLE HOLDER

Delicate version of the popular Crile-Wood needle holder. With TC jaws. Intended for 3/0 suture and finer.

| ART. NUMBER | PATTERN | LENGTH | TC |
|---------------|----------|----------|----|
| 104-15566-151 | Delicate | 15 cm 6" | • |
| 104-15566-181 | Delicate | 18 cm 7" | • |

OLSEN-HEGAR STILLE NEEDLE HOLDER

Popular needle holder with scissors for suturing and cutting with the same instrument. Intended for 0 suture and finer, 3+3 ratchets.

| ART. NUMBER | JAWS | LENGTH | TC |
|---------------|----------|---------------|----|
| 104-15590-121 | Scissors | 12 cm 4.75" | • |
| 104-15590-131 | Scissors | 13.5 cm 5.25" | • |



Art. Number: 104-15590-131



Art. Number: 104-15140



104-15140

RAPP-STILLE KNOT-HOLDING FORCEPS

A 13 cm knot-holding forceps with smooth jaws. Ideal instrument to use as a tool for holding the suture knot in narrow areas.

| ART. NUMBER | LENGTH |
|-------------|----------|
| 104-15140 | 13 cm 5" |

“I still have STILLE needle holders that date back from the 1970's and they are performing just as well as the new recent ones I purchased.”

- Dr. Mark L. Jewell, M.D.
Plastic Surgeon

“STILLE stands for impeccable quality, they're dependable and the durability is just outstanding.”

- Dr. Mark L. Jewell, M.D.
Plastic Surgeon

**For a complete presentation
of all STILLE needle holders,**

Scan the QR code to visit our website and
see our product catalogue



950-G104-0001-EN-01 © 2024-02-21 STILLE. All rights reserved.

Surgical perfection. For life.



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30 Torshälla, Sweden
☎ +46 8 588 580 00
✉ info@stille.se
🌐 www.stille.se

STILLE SURGICAL INC.

840 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614

DISTRIBUTED BY: