



Surgical perfection. For life.

STILLE Medstone Pain Management / ERCP / EP


MEDTECH OUTLOOK TOP 10
MEDICAL IMAGING
SOLUTION PROVIDERS - 2018

High Quality Hospital Grade Surgical Imaging Tables

The STILLE Medstone surgical imaging tables contribute to enhancing efficiency in the OR. Representing value priced solutions with quality and reliability.

Standard Configuration and Specifications

- Contoured tabletop 22-82" with facial cutout (Rectangular options available)
- Heavy duty steel base frame
- Four high quality hospital grade locking casters
- Mattress pad 2" thick with "Softcare™" technology
- Hand control for all table movements
- Patient safety straps (set of 3)
- Metal free imaging 58"
- Patient weight capacity 500lbs (weight upgrade options available)
- Parts and labor warranty 3/2 year
- 110VAC, 60Hz



Medstone5 PM with infectious control covers

The STILLE Medstone pain management series offers a broad range of motorized movements manufactured with high quality materials for durability and reliability.

1-5 movement options available

- Longitudinal
- Transverse
- Elevation
- Trendelenburg
- Lateral tilt



Optional arm board



Optional 3 movement foot control shown

STILLE Medstone Pain Tables

STILLE Medstone Table Model	Part #	# Movements	Elevation Range	Lateral Tilt	Trendelenburg	Longitudinal Travel	Infectious Control Covers	Transverse Travel
Medstone5 PM(H) TM5(H)	500-1001	5	26-40"	15°±	15°±	16"	Yes	3"±
Medstone4 PM(H) TM4(H)	500-1003	4	25-39"	15°±	15°±	16"	Yes	x
Medstone4 PM TM4ER	500-1004	4	29-43"	15°±	15°±	16"	x	x
Medstone3 PM TM3ER	500-1005	3	25-39"	x	15°±	16"	x	x
Medstone3 PM TM3MP	500-1006	3	26-40"	15°±	15°±	x	x	x
Medstone2 PM TM2ER	500-1009	2	23-37"	x	x	16"	x	x
Medstone2 PM TT2	500-1013	2	25-39"	x	15°±	x	x	x
Medstone2 PM TL2	500-1011	2	25-39"	15°±	x	x	x	x

One movement table also available: **Medstone1 PM, 500-1015** (Elevation range **22-36"**).



Height adjustment



Trendelenburg



Lateral tilt



Longitudinal

STILLE Medstone Pain Table Options

Part number	Description
505-1101	Rectangular top 22x82" WITHOUT facial cutout
505-1102	Rectangular top 22x82" WITH facial cutout
505-1104	Patient weight upgrade 750lbs w/rectangular top 26x82" WITHOUT facial cutout
505-1105	Patient weight upgrade 750lbs w/rectangular top 26x82" WITH facial cutout
505-1110	Dual patient arm board. Slides under mattress with two arm supports
505-1112	Combination patient arm board with pad and attachment clamp. One (1) quantity
505-1142	One (1) movement foot control for tables with elevation (fits table 500-1015)
505-1106	Two (2) movement foot control for tables with elevation and longitudinal travel (fits table 500-1009)
505-1140	Two (2) movement foot control for tables with elevation and lateral tilt (fits table 500-1011)
505-1141	Two (2) movement foot control for tables with elevation and Trendelenburg (fits table 500-1013)
505-1107	Three (3) movement foot control for tables with elevation, Trendelenburg, and lateral tilt (fits table 500-1006)
505-1138	Three (3) movement foot control for tables with elevation, Trendelenburg, longitudinal travel (fits table 500-1005)
505-1108	Four (4) movement foot control for tables with elevation, Trendelenburg, lateral tilt, longitudinal travel (fits table 500-1003 and 500-1004)
505-1109	Five (5) movement foot control for tables with elevation, Trendelenburg, lateral tilt, longitudinal travel, transverse travel (fits table 500-1001)
505-1137	Extension cable hand control
535-1770	Radiation shield kit, 6 single panels with Velcro attachment (grey only)

Specification tolerance: ±1°, ±.5"



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30
Torshälla, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
sales@stille.se

www.stille.se

STILLE Surgical Inc.

915 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614
sales.us@stille.se