



Surgical perfection. For life.

STILLE Medstone PM Electrophysiology/ Spinal/ Pain Management

High Quality Hospital Grade Surgical Imaging Tables 500lb/227kg weight capacity

The STILLE Medstone c-arm tables contribute to enhancing efficiency while representing an affordable solution with quality and reliability provided by Swedish technology.




Photo shows Medstone5 PM (500-1001)

Build the Table according to Your Needs

Choose from 2 to 5 of these available movements:

- Height adjustment
- Longitudinal travel
- Lateral travel
- Trendelenburg
- Lateral tilt

STILLE Medstone5 PM Surgical Imaging Table

5-Movement c-arm table for electrophysiology, spinal, pain management and any other applications that require good quality imaging, high performance, reliability and ease of use at affordable costs.

- Motorized movements
- Patient weight capacity 500lb (227kg)
- Carbon fiber table top with 1mm Aluminium Equivalence
- Auto horizontal and center home positioning
- Four hospital grade locking casters for great mobility
- Metal free imaging 58" (147cm)

Technical Specifications for a 5-Movement Standard Configuration

- Motorized height adjustment 26.78"-39.76" (65.5-101cm)
- $\pm 15^\circ$ motorized Trendelenburg tilt
- $\pm 15^\circ$ motorized lateral tilt
- 15 3/4" (40cm) motorized longitudinal travel
- 6" (15cm) motorized lateral travel with auto centering
- Patient weight capacity 500lb (227kg)
- 230Vac, 3,4A, 50 Hz and 120Vac, 6.5A, 60 Hz with 230Vac or 120Vac style line cord
- Infectious control covers
- 2" (5cm) table mattress and 3 patient straps
- Wired hand control

Table highlights

STILLE Medstone5 PM (5-movement) table top dimensions 22"x 82" (56x208cm)

Item/ Movement	Description
500-1001	STILLE Medstone5 PM surgical imaging table with carbon fiber table top with facial cutout
Motorized Elevation Travel	minimum height 26.78" (65.5cm) to maximum height 39.76" (101cm)
Motorized Trendelenburg Tilt	± 15°
Motorized Lateral Tilt	± 15°
Motorized Longitudinal Travel	15 ¾" (40cm)
Motorized Lateral Travel	6" (15cm)

STILLE Medstone4 PM (4-movement) table top dimensions 22"x82" (56x208cm)

Item/ Movement	Description
500-1004	STILLE Medstone4 PM surgical imaging table with carbon fiber table top with facial cutout
Motorized Elevation Travel	minimum height 28.74" (73cm) to maximum height 42.52" (108cm)
Motorized Trendelenburg Tilt	± 15°
Motorized Longitudinal Travel	15 ¾" (40cm)
Motorized Lateral Tilt	± 15°

STILLE Medstone3 PM (3-movement) table top dimensions 22"x82" (56x208cm)

Item/ Movement	Description
500-1006	STILLE Medstone3 PM surgical imaging table with carbon fiber table top with facial cutout
Motorized Elevation Travel	minimum height 25.39" (64.5cm) to maximum height 40.16" (102cm)
Motorized Trendelenburg Tilt	± 15°
Motorized Lateral Tilt	± 15°

Various Table Tops available



A=PM table standard top (contoured with facial cutout-for prone patient positioning) with optional arm board 505-1112.



Alternative table top (rectangular without facial cutout-for supine patient positioning) available in standard and wide version; for part numbers see page three.

Customize the table movements according to your needs

Model	Medstone 5 PM	Medstone 4 PM	Medstone 4 PM	Medstone 3 PM	Medstone 3 PM	Medstone 2 PM	Medstone 2 PM	Medstone 2 PM
Part number	500-1001	500-1003	500-1004	500-1005	500-1006	500-1009	500-1011	500-1013
Height adjustment	26.78"-39.76" (65.5-101cm)	25.39"-39.37" (64.5-100cm)	28.74"-42.52" (73-108cm)	25"-38.78" (63.5-98.5cm)	25.39"-40.16" (64.5-102cm)	23.03"-37.00" (58.5-94cm)	25.39"-39.37" (64.5-100cm)	22.44"-36.22" (57-92cm)
Trendelenburg	± 15°	± 15°	± 15°	± 15°	± 15°	x	x	± 15°
Lateral tilt	± 15°	± 15°	± 15°	x	± 15°	x	± 15°	x
Longitud. travel	40cm	40cm	40cm	40cm	x	40cm	x	x
Lateral travel	15cm	x	x	x	x	x	x	x
Auto horizontal & center home positioning	✓	✓	x	x	x	x	x	x
Infectious control covers	✓	✓	x	x	x	x	x	x
Table top	B	B	A	A	A	A	A	A

Specification tolerance: ±1°, ±.5" (1.27cm)

A=Contoured carbon fibre table top 22"x82" (56x208cm) with facial cutout. This table top shape is recommended for pain management.

B= Same as A plus a 11.97" (30.4cm) long side rail fixed to the table top (as shown in the photo on the first page of this brochure).

C=Optional rectangular carbon fibre table tops (available in different sizes).



Height adjustment



Trendelenburg



Lateral tilt



Longitudinal travel

STILLE Medstone optional accessories

Part number	Description
505-1101	Rectangular top 22"x82" (56x208cm) WITHOUT facial cutout*
505-1102	Rectangular top 22"x82" (56x208cm) WITH facial cutout
505-1104	Wide rectangular top 26"x82" (66x208cm) WITHOUT facial cutout and with upgraded weight capacity 750lbs (340kg)*
505-1105	Wide rectangular top 26"x82" (66x208cm) WITH facial cutout and upgraded weight capacity 750lbs (340kg)
505-1179	2 detachable side rails length 31.69" (80.5cm) for all models with tops 22"x82" (56x208cm) except tables with table top B
505-1175	2 detachable side rails length 44.29" (112.5cm) only for the models 500-1001 and 500-1003 with table top 22"x82" (56x208cm)
505-1181	2 detachable side rails length 31.69" (80.5cm) for all tables with wide table tops 26"x82" (66x208cm)
514-50-3SR025	Universal adjustable arm board with clamp (side rail mounted)
505-1112	Arm board with pad and clamp for attachment directly to the table top
505-1110	Dual arm board with pads; to be placed across the table top; no clamps or side rails needed
514-100401	IV-pole with clamp (side rail mounted)
514-100311	Anesthesia frame with clamps (side rail mounted)

* Rectangular table tops without facial cutouts are recommended for electrophysiology and spinal.

STILLE Medstone optional accessories

Part number	Description
535-1770	X-ray protective kit grey, 6 single panels; attachment to table top with self-adhesive Velcro © strips
535-1771	X-ray protective kit blue, 6 single panels; attachment to table top with self-adhesive Velcro © strips
535-1772	X-ray protective loop grey, 5 connected panels, attachment to side rails or table top with self-adhesive Velcro © strips
535-1773	X-ray protective loop blue, 5 connected panels, attachment to side rails or table top with self-adhesive Velcro © strips
505-1106	Foot control for Medstone2 PM (500-1009)
505-1140	Foot control for Medstone2 PM (500-1011)
505-1141	Foot control for Medstone2 PM (500-1013)
505-1107	Foot control for Medstone3 PM (500-1006)
505-1138	Foot control for Medstone3 PM (500-1005)
505-1108	Foot control for Medstone4 PM (500-1003 and 500-1004)
505-1109	Foot control for Medstone5 PM (500-1001)

STILLE – Surgical perfection since 1841

STILLE is one of the oldest medical device companies in the world founded at Karolinska University Hospital in Sweden back in 1841. At STILLE, we develop, produce and sell high-quality surgical instruments and surgical imaging tables to physicians worldwide.

STILLE's offering consists of a strong brand and products renowned for their superior quality. We manage our business in two distinct business areas—Surgical Instruments and Surgical imaging tables.

Surgical instruments

In our Surgical Instruments business area, we have a complete offering of handcrafted, forged instrumentation for all fields of open surgery, with a special emphasis on cardiovascular surgery, and reconstructive and aesthetic plastic surgery. Our instruments feature delicate design, industry-leading durability and a unique feel, which has made STILLE a well-known and highly respected supplier by surgeons throughout the world.

Surgical imaging tables

In our Surgical imaging tables business area, we offer specialty tables for minimally invasive procedures, with pricing, safety and features optimized for both in- and outpatient settings. The STILLE imagiQ2 surgical imaging table has the highest translucency in the world, and is optimized for advanced mobile Hybrid OR settings as well as peripheral and interventional procedures, enhancing crystal imaging while reducing radiation exposure as well as increasing efficiency in the OR.

The STILLE Medstone line of surgical imaging tables range from floating to fully motorized tables for vascular, electrophysiology, gastro (ERCP), pain management and urology to suit your needs, with reliability and optimal performance at an affordable cost at the highest standards of Swedish technology.



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30
Torshälla, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
sales@stille.se

www.stille.se

STILLE Surgical Inc.

915 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614
sales.us@stille.se