




Surgical perfection. For life.

# Year codes STILLE Instruments

All STILLE instruments are marked with a year code. The code for respective year is found in the table.



On the STILLE scissors the year code is inscribed on the inside of one of the shanks.

## Year Codes STILLE instruments

YEAR	CODE
1938	8
1939	9
1940	0
1941	1
1942	2
1943-1948	Unstamped
1949	•
1950	•
1951	1
1952	1
1953	Unstamped
1954	Unstamped
1955	E9
1956	F9
1957	G9
1958	H9
1959	I9
1960	K9
1961	L9

YEAR	CODE
1962	M9
1963	N9
1964	O9
1965	P9
1966	Q9
1967	R9
1968	S9
1969	T9
1970	U9
1971	V9
1972	X9
1973	Y9
1974	Z9
1975	A10
1976	B10
1977	C10
1978	D10
1979	E10
1980	F10

YEAR	CODE
1981	G10
1982	H10
1983	A
1984	B
1985	C
1986	D
1987	E
1988	F
1989	G
1990	H
1991	I
1992	K
1993	L
1994	M
1995	O
1996	P
1997	R
1998	S
1999	T

YEAR	CODE
2000	U
2001	X
2002	Y
2003	Z
2004	4
2005	5
2006	6
2007	7
2008	8
2009	9
2010	10
2011	11
2012	12
2013	13
2014	14
2015	15
2016	16